

sigvaris

SIGVARIS
GROUP



Healthcare Professional
Guide of Prescription
Products

Healthcare Professional Guide of Prescription Products

Segment	Product Name
Circular knit	Essential COMFORTABLE
	Essential THERMOREGULATING
	Active WORK WEAR
	Active LEISURE
	TRADITIONAL
	ULCER X Kit
	ULCER X Liner
	Essential SEMITRANSSPARENT
	Style OPAQUE
	Style SEMITRANSSPARENT
	Style TRANSPARENT
	Style PATTERNS
	Essential COTON
	Essential MICROFIBRE
	Style COLOURS
	Active MASCULINE
Inelastic compression wraps	COMPREFLEX Standard - Calf
	COOLFLEX Standard - Calf
	COMPREFLEX Standard - Thigh
	COMPREFLEX Standard - Knee
	COMPREBOOT Standard
Upper limb	COMPREFLEX Standard - Arm
	ADVANCE Armsleeve
	TRADITIONAL Armsleeve
Donning and doffing aids	DOFF N'DONNER ROLLY
	DOFF N'DONNER CONE
	MELANY

sigvaris

We're here to help

Customer Services

Need help? Just ask - we're here to make it easy.

- Our friendly product experts have over 150 years of combined experience in compression therapy - no question is too small
- Made to measure support: we'll help you get the perfect fit for your patients
- We're here Monday to Thursday from 9am to 5pm and Friday from 9am to 4pm - just a call or email away

For any queries, please refer to the customer service contact details at the bottom of each page.



Call us on
01264 326 666
or email us on
**customerservices@
sigvaris.com**

Ensuring the correct fit will allow the legs to receive the full therapeutic benefits of compression.

That's why Sigvaris provides guidance on the best measuring points to achieve an optimum fit.

General guidelines for measuring

- Measurements should ideally be taken as early in the morning as possible, when swelling is minimal.
- Where possible, measure against bare skin to ensure accuracy.
- The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images to the right for examples).
- Circumference measurements should be taken by pulling the tape around the leg so that it is snug but not constricting.
- Length measurements should be taken from the ground up along a straight tape measure up the outside of the leg.




The diagrams below illustrate the available measurements and indicate where measurements should be taken based on the selected product style.

Circular knit measuring points

Definition potential of the measuring points for circular knit compression

cB	Measure the circumference of the ankle at the narrowest point
cC	Measure the circumference of the calf at the fullest point
cG	Measure the circumference of the widest part of the thigh
cH	Measure the circumference at the hip
cT	For tights (panty), measure the circumference at the waist
ℓD	Measure the full calf length from the floor to 2cm below the right-angled bend of the knee
ℓG	Measure the leg length from the floor to 2cm below the gluteal fold


c - Circumference ℓ - Length



Inelastic compression wrap measuring points

Definition potential of the measuring points for inelastic compression wraps

B	Measure the ankle circumference just above the narrowest part of the ankle
B1	Measure the circumference just below the calf
C	Measure the calf circumference at the fullest part
D	Measure 5cm below the top of the fibula
E	Measure 5cm above the top of the patella
F	Measure the circumference mid way up the thigh
G	Measure the circumference 2cm below the gluteal fold
L1	Determine the calf length by measuring from the ankle bone to the bend of the knee
L2	Determine the thigh length by measuring from the top of the patella to the gluteal fold
Y	Measure around the ankle bend and heel
A	Measure around the metatarsal heads
K	Measuring from the inside of the foot from the heel to the big toe metatarsal



Scan the QR code to view our 'How to' measure videos

For all **Made to Measure (MTM)** and **Upper limb** measuring points, please use the QR code to browse the literature on our website.



Scan the QR code to view all product literature

Access to compression therapy should be straightforward and efficient to support the timely management of venous and lymphatic conditions. In the UK, there are several established routes for ordering medical compression garments.

Each route is designed to meet different clinical and logistical needs, helping ensure patients receive the appropriate product.

Using a Dispensing Appliance Contractors (DAC)

DACs are specialised services that provide medical appliances, including compression garments, directly to patients.

They work with prescriptions from GPs and offer home delivery services. Examples of DACs include **Daylong** and **Patient Choice**.

Their aim is to offer a prompt service, home delivery, patient support and reduce appliance wastage.



NHS Supply chain

Many Sigvaris products are available through NHS Supply Chain, making them easily orderable by hospitals and community health services across the UK.

To order through this route, first identify the product and size you need. Then, use the corresponding 5-7 digit code listed in the Product Code catalogue to place your order.




Scan the QR code to view our product code catalogue

Introducing Sigvaris Select

Your new tool to find Sigvaris garments with precision, speed and confidence.

Why guess when you can Select - try it now



Need help finding a product on EMIS, SystemOne or with GP support?

Call us on 01264 326 666 or email us on customerservices@sigvaris.com



Scan the QR code to read our Ordering compression garments and accessories guide

Pharmacy

Patients can take their prescription from their prescriber directly to a local pharmacy, where the pharmacy can order Sigvaris products. The prescription must include:

- Company (e.g. Sigvaris)
- Product name
- mmHg (Compression class e.g. RAL class 1, 2, or 3)
- Style (product style e.g. below knee, thigh, tights etc.)
- Size (size of the garment e.g. small, medium, large etc.)
- Toe options (if required)
- Colour
- Number of pairs needed

GP surgeries

Healthcare professionals can initiate the prescription directly through EMIS or similar systems at the GP surgery. Sigvaris supports EMIS search tips for efficient prescribing. Clinicians can also access Sigvaris product codes to ensure accurate garment selection and ordering.



Searching for prescription products using EMIS

Below are some tips for using the EMIS online system to narrow your list of choices.

- Searching for Sigvaris products occurs in the 'Add Drug' screen within a patient record.
- The search takes place by typing into the name box - firstly with the EMIS tip product name.
- The EMIS web based system search function works best when you type in abbreviated versions of the words.
- These abbreviated words do not need to be typed in any particular order.
- Use the keyboard down arrow to scroll through the list. This will allow you to see the full product description.

Sigvaris compression - in a class of its own

SIGVARIS medical compression stockings are available in AFNOR and RAL compression standards with stringent defined procedures on how to produce and test.

AFNOR

Products generally have a soft feel and offer low to medium levels of medical compression with the range anatomically shaped for men and women.

RAL

Products can offer higher compression levels and very often are unisex garments with some products consisting of firmer, more supportive materials.

		Light					Strong
Product		Compression classes	10 mmHg	20 mmHg	30 mmHg	40 mmHg	50 mmHg
Circular knit	Essential COMFORTABLE	RAL 1-2		18-21	23-32 mmHg		
	Essential THERMOREGULATING	RAL 1-3		18-21	23-32 mmHg	34-46 mmHg	
	Active WORK WEAR	AFNOR 2	15-20 mmHg				
	Active LEISURE	AFNOR 2	15-20 mmHg				
	TRADITIONAL	RAL 2-3			23-32 mmHg	34-46 mmHg	
	ULCER X Kit	>45 mmHg					>45 mmHg
	ULCER X Liner	15-20 mmHg		15-20 mmHg			
	Essential SEMITRANSSPARENT	AFNOR 1-3	10-15 mmHg	15-20 mmHg	20-36 mmHg		
	Style OPAQUE	AFNOR 2		15-20 mmHg			
	Style SEMITRANSSPARENT	RAL 1-2		18-21	23-32 mmHg		
	Style TRANSPARENT	AFNOR 2		15-20 mmHg			
	Style PATTERNS	AFNOR 2		15-20 mmHg			
	Essential COTON	AFNOR 2-3		15-20 mmHg	20-36 mmHg		
	Essential MICROFIBRE	AFNOR 2		15-20 mmHg			
Style COLOURS	AFNOR 2		15-20 mmHg				
Active MASCULINE	RAL 2			23-32 mmHg			
Made to Measure (MTM)	MTM Essential COMFORTABLE	RAL 1-2		18-21	23-32 mmHg		
	MTM Essential THERMOREGULATING	RAL 1-3		18-21	23-32 mmHg	34-46 mmHg	
	MTM TRADITIONAL	RAL 2-4			23-32 mmHg	34-46 mmHg	>49 mmHg
	MTM Style SEMITRANSSPARENT	RAL 1-2		18-21	23-32 mmHg		
	MTM OPTIFORM Flex - Flat knit	RAL 1-2		18-21	23-32 mmHg		
	MTM OPTIFORM Hold - Flat knit	RAL 1-3		18-21	23-32 mmHg	34-46 mmHg	
Inelastic compression wraps	COMPREFLEX Standard - Calf	20-50 mmHg			20-50 mmHg Variable		
	COMPREFLEX Standard - Thigh	20-50 mmHg			20-50 mmHg Variable		
	COMPREFLEX Standard - Knee	20-50 mmHg			20-50 mmHg Variable		
	COMPREBOOT Standard	30-40 mmHg			30-40 mmHg Variable		
	COOLFLEX Standard - Calf	20-50 mmHg			20-50 mmHg Variable		
Upper limb	COMPREFLEX Standard - Arm	20-30 mmHg		20-30 mmHg Variable			
	ADVANCE Armsleeve	14-18 mmHg 20-25 mmHg 30-40 mmHg	14-18	20-25 mmHg	30-40 mmHg		
	TRADITIONAL Armsleeve	RAL 2			23-32 mmHg		

Let's Measure Together

We have CPD accredited workshops available to support making measuring simple for Sigvaris compression garments.

Sessions accredited by

Scan the QR code to book a session



Circular knit

Products available

Essential COMFORTABLE	14-15
Essential THERMOREGULATING	16-17
Active WORK WEAR	18-19
Active LEISURE	20-21
TRADITIONAL	22-23
ULCER X Kit	24-25
ULCER X Liner	26-27
Essential SEMITRANSSPARENT	28-29
Style OPAQUE	30-31
Style SEMITRANSSPARENT	32-33
Style TRANSPARENT	34-35
Style PATTERNS	36-37
Essential COTON	38-39
Essential MICROFIBRE	40-41
Style COLOURS	42-43
Active MASCULINE	44-45

Available in MTM

Offers tailored sizing, colour options, and style choices, ensuring the right product with the right fit.

Find out more in the Made to Measure section

Material composition

All styles

65% Polyamide, 35% Elastane (content can vary slightly by compression choice)
Latex free



Scan here for this product's measuring form

Essential COMFORTABLE

The modern fabric provides a soft and smooth feel against the skin while offering breathability and moisture-wicking to aid skin care.

Compression class

18-21 mmHg RAL1
23-32 mmHg RAL2

Product benefits

Soft, Easy donning & doffing, Moisture-wicking, Available in MTM

Suitable for

- Venous insufficiency
- Early stages of lymphoedema / chronic oedema
- Active leg ulcers
- Mixed disease
- Lipoedema (lipalgia)
- Prevention of leg ulcers and oedema



Knobbed grip top band



Thigh with grip top in Black



Thigh with waist attachment in Savannah in Open toe



Tights (panty) in Caramel

Below knee in Caramel with Open toe

EMIS search tips

Product name
Essential COMFORTABLE

SIG ESS COMF
EMIS TIP

Compression class (mmHg level)
Class 1 18 Class 2 23
EMIS TIP EMIS TIP

Sizing chart


All measurements are shown in cm

Below knee, thigh and tights (panty)

	XXS		XS		S		M		L		XL	
Ankle (cB)	18-20		20-22		22-24		24-26		26-29		29-32	
	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS
Calf (cC)	25-32	29-36	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49	42-49	43-51
Thigh (cG) with grip top	42-52	47-56	47-56	52-61	52-61	56-65	56-65	61-71	61-71	66-76	66-76	73-85
Thigh (cG) (For thigh w/ attachment and tights only)	42-54	47-59	47-59	52-64	52-64	56-68	56-68	61-74	61-74	66-79	66-79	76-85
Hip (cH) (For thigh w/ attachment and tights only)	75-90	80-95	80-95	85-100	85-100	95-110	95-110	105-120	105-120	115-130	115-130	125-140
Waist (cT) (For thigh w/ attachment and tights only)	60-75	65-80	65-80	75-90	75-90	85-100	85-100	95-110	95-110	105-120	105-120	115-130
EMIS TIP	SMA	SMA PLU	SMA	SMA PLU	SMA	SMA PLU	MED	MED PLU	LAR	LAR PLU	LAR	LAR PLU

Length

	Short	Normal	Long
Below knee (lD)	-	35-40	41-45
Thigh (lG)	56-64	65-74	75-84
Tights (panty) (lG)	-	65-74	75-84
EMIS TIP	SHO	NOR	LON



Styles and colours

●●● Closed toe CLO
○○○ Open toe OPE

		Class 1			Class 2		
		Black	Caramel	Savannah	Black	Caramel	Savannah
Below knee	BEL	●●	●●	●●	●●	●●	●●
Thigh without grip top	THI	●●	●●	●●	●●	●●	●●
Thigh with grip top	THI GRI	●●	●●	●●	●●	●●	●●
Thigh with waist attachment (Left and right)	THI WAI LEF	-	-	-	●●	●●	●●
	THI WAI RIG	-	-	-	●●	●●	●●
Tights (panty)	PAN	●●	●●	●●	●●	●●	●●
Tights (panty) extra	PAN	●●	●●	●●	●●	●●	●●
Tights (panty) body form	BOD	-	-	-	●●	●●	●●
EMIS TIP		BLK	CAR	SAV	BLK	CAR	SAV



Scan here for this product's ordering codes

Essential COMFORTABLE



Product codes

Below knee


FOR THIS PRODUCT - Please add '-CM' to the end of your product codes when ordering

XXS	STD	Normal	Class 1						Class 2						
			Black		Caramel		Savannah		Black		Caramel		Savannah		
			●	○	●	○	●	○	●	○	●	○	●	○	
		Long	68106	68098	68074	68066	68138	68114	69522	69514	69490	69482	69554	69530	
		PLUS	Normal	68109	68101	68077	68069	68149	68125	69525	69517	69493	69485	69565	69541
		Long	68110	68102	68078	68070	68150	68126	69526	69518	69494	69486	69566	69542	
			STD	Normal	47513	47417	47481	47385	68139	68115	39144	39128	39080	39064	69555
Long	47514			47418	47482	47386	68140	68116	39145	39129	39081	39065	69556	69532	
		PLUS	Normal	47521	47425	47489	47393	68151	68127	39152	39136	39088	39072	69567	69543
			Long	47522	47426	47490	47394	68152	68128	39153	39137	39089	39073	69568	69544
		STD	Normal	47515	47419	47483	47387	68141	68117	39146	39130	39082	39066	69557	69533
			Long	47516	47420	47484	47388	68142	68118	39147	39131	39083	39067	69558	69534
		PLUS	Normal	47523	47427	47491	47395	68153	68129	39154	39138	39090	39074	69569	69545
			Long	47524	47428	47492	47396	68154	68130	39155	39139	39091	39075	69570	69546
		STD	Normal	47517	47421	47485	47389	68143	68119	39148	39132	39084	39068	69559	69535
			Long	47518	47422	47486	47390	68144	68120	39149	39133	39085	39069	69560	69536
		PLUS	Normal	47525	47429	47493	47397	68155	68131	39156	39140	39092	39076	69571	69547
			Long	47526	47430	47494	47398	68156	68132	39157	39141	39093	39077	69572	69548
		STD	Normal	47519	47423	47487	47391	68145	68121	39150	39134	39086	39070	69561	69537
			Long	47520	47424	47488	47392	68146	68122	39151	39135	39087	39071	69562	69538
		PLUS	Normal	47527	47431	47495	47399	68157	68133	39158	39142	39094	39078	69573	69549
			Long	47528	47432	47496	47400	68158	68134	39159	39143	39095	39079	69574	69550
		STD	Normal	68107	68099	68075	68067	68147	68123	69523	69515	69491	69483	69563	69539
			Long	68108	68100	68076	68068	68148	68124	69524	69516	69492	69484	69564	69540
		PLUS	Normal	68111	68103	68079	68071	68159	68135	69527	69519	69495	69487	69575	69551
			Long	68112	68104	68080	68072	68160	68136	69528	69520	69496	69488	69576	69552

●●● Closed toe ○○○ Open toe

Product codes

Thigh without grip top




FOR THIS PRODUCT - Please add '-CM' to the end of your product codes when ordering

XXS	STD	Short	Class 1						Class 2						
			Black		Caramel		Savannah		Black		Caramel		Savannah		
			●	○	●	○	●	○	●	○	●	○	●	○	
		Normal	68197	68185	68173	68161	68245	68209	69757	69745	69733	69721	69805	69769	
		Long	68199	68187	68175	68163	68247	68211	69759	69747	69735	69723	69807	69771	
		PLUS	Short	68203	68191	68179	68167	68263	68227	69763	69751	69739	69727	69823	69787
			Normal	68204	68192	68180	68168	68264	68228	69764	69752	69740	69728	69824	69788
		Long	Short	68205	68193	68181	68169	68265	68229	69765	69753	69741	69729	69825	69789
			Normal	68206	68194	68182	68170	68266	68230	69766	69754	69742	69730	69826	69790
		STD	Short	47793	47649	47745	47601	68248	68212	39424	39400	39328	39304	69808	69772
			Normal	47794	47650	47746	47602	68249	68213	39425	39401	39329	39305	69809	69773
		Long	Short	47795	47651	47747	47603	68250	68214	39426	39402	39330	39306	69810	69774
			Normal	47796	47652	47748	47604	68251	68215	39427	39403	39331	39307	69811	69775
		PLUS	Short	47805	47661	47757	47613	68266	68230	39436	39412	39340	39316	69826	69790
			Normal	47806	47662	47758	47614	68267	68231	39437	39413	39341	39317	69827	69791
		Long	Short	47807	47663	47759	47615	68268	68232	39438	39414	39342	39318	69828	69792
			Normal	47808	47664	47760	47616	68269	68233	39439	39415	39343	39319	69829	69793
		STD	Short	47796	47652	47748	47604	68251	68215	39427	39403	39331	39307	69811	69775
			Normal	47797	47653	47749	47605	68252	68216	39428	39404	39332	39308	69812	69776
		Long	Short	47798	47654	47750	47606	68253	68217	39429	39405	39333	39309	69813	69777
			Normal	47799	47655	47751	47607	68254	68218	39430	39406	39334	39310	69814	69778
		PLUS	Short	47808	47664	47760	47616	68269	68233	39439	39415	39343	39319	69829	69793
			Normal	47809	47665	47761	47617	68270	68234	39440	39416	39344	39320	69830	69794
		Long	Short	47810	47666	47762	47618	68271	68235	39441	39417	39345	39321	69831	69795
			Normal	47811	47667	47763	47619	68272	68236	39442	39418	39346	39322	69832	69796
		STD	Short	47799	47655	47751	47607	68254	68218	39430	39406	39334	39310	69814	69778
			Normal	47800	47656	47752	47608	68255	68219	39431	39407	39335	39311	69815	69779
		Long	Short	47801	47657	47753	47609	68256	68220	39432	39408	39336	39312	69816	69780
			Normal	47812	47668	47764	47620	68273	68237	39443	39419	39347	39323	69833	69797
		PLUS	Short	47811	47667	47763	47619	68272	68236	39442	39418	39346	39322	69832	69796
			Normal	47812	47668	47764	47620	68273	68237	39443	39419	39347	39323	69833	69797
		Long	Short	47813	47669	47765	47621	68274	68238	39444	39420	39348	39324	69834	69798
			Normal	47814	47670	47766	47622	68275	68239	39445	39421	39349	39325	69835	69799
		STD	Short	47802	47658	47754	47610	68257	68221	39433	39409	39337	39313	69817	69781
			Normal	47803	47659	47755	47611	68258	68222	39434	39410	39338	39314	69818	69782
		Long	Short	47804	47660	47756	47612	68259	68223	39435	39411	39339	39315	69819	69783
			Normal	47815	47671	47767	47623	68276	68240	39446	39422	39350	39326	69836	69800
		PLUS	Short	47814	47670	47766	47622	68275	68239	39445	39421	39349	39325	69835	69799
			Normal	47815	47671	47767	47623	68276	68240	39446	39422	39350	39326	69836	69800
		Long	Short	47816	47672	47768	47624	68277	68241	39447	39423	39351	39327	69837	69801
			Normal	68200	68188	68176	68164	68260	68224	69760	69748	69736	69724	69820	69784
		STD	Normal	68201	68189	68177	68165	68261	68225	69761	69749	69737	69725	69821	69785
			Long	68202	68190	68178	68166	68262	68226	69762	69750	69738	69726	69822	69786
		PLUS	Short	68206	68194	68182	68170	68278	68242	69766	69754	69742	69730	69838	69802
			Normal	68207	68195	68183	68171	68279	68243	69767	69755	69743	69731	69839	69803
		Long	Short	68208	68196	68184	68172	68280	68244	69768	69756	69744	69732	69840	69804

●●● Closed toe ○○○ Open toe

Product codes

Thigh with grip top




FOR THIS PRODUCT - Please add '-CM' to the end of your product codes when ordering

			Class 1						Class 2						
			Black		Caramel		Savannah		Black		Caramel		Savannah		
			●	○	●	○	●	○	●	○	●	○	●	○	
XXS	STD	Short	68605	68569	68389	68353	68533	68497	70165	70129	69949	69913	70093	70057	
		Normal	68606	68570	68390	68354	68534	68498	70166	70130	69950	69914	70094	70058	
		Long	68607	68571	68391	68355	68535	68499	70167	70131	69951	69915	70095	70059	
	PLUS	Short	68623	68587	68407	68371	68551	68515	70183	70147	69967	69931	70111	70075	
		Normal	68624	68588	68408	68372	68552	68516	70184	70148	69968	69932	70112	70076	
		Long	68625	68589	68409	68373	68553	68517	70185	70149	69969	69933	70113	70077	
	XS	STD	Short	68608	68572	68392	68356	68536	68500	70168	70132	69952	69916	70096	70060
			Normal	68609	68573	68393	68357	68537	68501	70169	70133	69953	69917	70097	70061
			Long	68610	68574	68394	68358	68538	68502	70170	70134	69954	69918	70098	70062
PLUS		Short	68626	68590	68410	68374	68554	68518	70186	70150	69970	69934	70114	70078	
		Normal	68627	68591	68411	68375	68555	68519	70187	70151	69971	69935	70115	70079	
		Long	68628	68592	68412	68376	68556	68520	70188	70152	69972	69936	70116	70080	
S	STD	Short	68611	68575	68395	68359	68539	68503	70171	70135	69955	69919	70099	70063	
		Normal	68612	68576	68396	68360	68540	68504	70172	70136	69956	69920	70100	70064	
		Long	68613	68577	68397	68361	68541	68505	70173	70137	69957	69921	70101	70065	
	PLUS	Short	68629	68593	68413	68377	68557	68521	70189	70153	69973	69937	70117	70081	
		Normal	68630	68594	68414	68378	68558	68522	70190	70154	69974	69938	70118	70082	
		Long	68631	68595	68415	68379	68559	68523	70191	70155	69975	69939	70119	70083	
M	STD	Short	68614	68578	68398	68362	68542	68506	70174	70138	69958	69922	70102	70066	
		Normal	68615	68579	68399	68363	68543	68507	70175	70139	69959	69923	70103	70067	
		Long	68616	68580	68400	68364	68544	68508	70176	70140	69960	69924	70104	70068	
	PLUS	Short	68632	68596	68416	68380	68560	68524	70192	70156	69976	69940	70120	70084	
		Normal	68633	68597	68417	68381	68561	68525	70193	70157	69977	69941	70121	70085	
		Long	68634	68598	68418	68382	68562	68526	70194	70158	69978	69942	70122	70086	
L	STD	Short	68617	68581	68401	68365	68545	68509	70177	70141	69961	69925	70105	70069	
		Normal	68618	68582	68402	68366	68546	68510	70178	70142	69962	69926	70106	70070	
		Long	68619	68583	68403	68367	68547	68511	70179	70143	69963	69927	70107	70071	
	PLUS	Short	68635	68599	68419	68383	68563	68527	70195	70159	69979	69943	70123	70087	
		Normal	68636	68600	68420	68384	68564	68528	70196	70160	69980	69944	70124	70088	
		Long	68637	68601	68421	68385	68565	68529	70197	70161	69981	69945	70125	70089	
XL	STD	Short	68620	68584	68404	68368	68548	68512	70180	70144	69964	69928	70108	70072	
		Normal	68621	68585	68405	68369	68549	68513	70181	70145	69965	69929	70109	70073	
		Long	68622	68586	68406	68370	68550	68514	70182	70146	69966	69930	70110	70074	
	PLUS	Short	68638	68602	68422	68386	68566	68530	70198	70162	69982	69946	70126	70090	
		Normal	68639	68603	68423	68387	68567	68531	70199	70163	69983	69947	70127	70091	
		Long	68640	68604	68424	68388	68568	68532	70200	70164	69984	69948	70128	70092	

●●● Closed toe ○○○ Open toe

Product codes

Thigh with waist attachment - Left




FOR THIS PRODUCT - Please add '-CM' to the end of your product codes when ordering

			Class 2						
			Black		Caramel		Savannah		
			●	○	●	○	●	○	
XXS	STD	Short	80186	79970	80114	79898	80150	79934	
		Normal	80187	79971	80115	79899	80151	79935	
		Long	80188	79972	80116	79900	80152	79936	
	PLUS	Short	80204	79988	80132	79916	80168	79952	
		Normal	80205	79989	80133	79917	80169	79953	
		Long	80206	79990	80134	79918	80170	79954	
	XS	STD	Short	80189	79973	80117	79901	80153	79937
			Normal	80190	79974	80118	79902	80154	79938
			Long	80191	79975	80119	79903	80155	79939
PLUS		Short	80207	79991	80135	79919	80171	79955	
		Normal	80208	79992	80136	79920	80172	79956	
		Long	80209	79993	80137	79921	80173	79957	
S	STD	Short	80192	79976	80120	79904	80156	79940	
		Normal	80193	79977	80121	79905	80157	79941	
		Long	80194	79978	80122	79906	80158	79942	
	PLUS	Short	80210	79994	80138	79922	80174	79958	
		Normal	80211	79995	80139	79923	80175	79959	
		Long	80212	79996	80140	79924	80176	79960	
M	STD	Short	80195	79979	80123	79907	80159	79943	
		Normal	80196	79980	80124	79908	80160	79944	
		Long	80197	79981	80125	79909	80161	79945	
	PLUS	Short	80213	79997	80141	79925	80177	79961	
		Normal	80214	79998	80142	79926	80178	79962	
		Long	80215	79999	80143	79927	80179	79963	
L	STD	Short	80198	79982	80126	79910	80162	79946	
		Normal	80199	79983	80127	79911	80163	79947	
		Long	80200	79984	80128	79912	80164	79948	
	PLUS	Short	80216	80000	80144	79928	80180	79964	
		Normal	80217	80001	80145	79929	80181	79965	
		Long	80218	80002	80146	79930	80182	79966	
XL	STD	Short	80201	79985	80129	79913	80165	79949	
		Normal	80202	79986	80130	79914	80166	79950	
		Long	80203	79987	80131	79915	80167	79951	
	PLUS	Short	80219	80003	80147	79931	80183	79967	
		Normal	80220	80004	80148	79932	80184	79968	
		Long	80221	80005	80149	79933	80185	79969	

●●● Closed toe ○○○ Open toe

Product codes

Thigh with waist attachment - Right




FOR THIS PRODUCT - Please add '-CM' to the end of your product codes when ordering

			Class 2					
			Black		Caramel		Savannah	
			●	○	●	○	●	○
XXS	STD	Short	80294	80078	80222	80006	80258	80042
		Normal	80295	80079	80223	80007	80259	80043
		Long	80296	80080	80224	80008	80260	80044
	PLUS	Short	80312	80096	80240	80024	80276	80060
		Normal	80313	80097	80241	80025	80277	80061
		Long	80314	80098	80242	80026	80278	80062
XS	STD	Short	80297	80081	80225	80009	80261	80045
		Normal	80298	80082	80226	80010	80262	80046
		Long	80299	80083	80227	80011	80263	80047
	PLUS	Short	80315	80099	80243	80027	80279	80063
		Normal	80316	80100	80244	80028	80280	80064
		Long	80317	80101	80245	80029	80281	80065
S	STD	Short	80300	80084	80228	80012	80264	80048
		Normal	80301	80085	80229	80013	80265	80049
		Long	80302	80086	80230	80014	80266	80050
	PLUS	Short	80318	80102	80246	80030	80282	80066
		Normal	80319	80103	80247	80031	80283	80067
		Long	80320	80104	80248	80032	80284	80068
M	STD	Short	80303	80087	80231	80015	80267	80051
		Normal	80304	80088	80232	80016	80268	80052
		Long	80305	80089	80233	80017	80269	80053
	PLUS	Short	80321	80105	80249	80033	80285	80069
		Normal	80322	80106	80250	80034	80286	80070
		Long	80323	80107	80251	80035	80287	80071
L	STD	Short	80306	80090	80234	80018	80270	80054
		Normal	80307	80091	80235	80019	80271	80055
		Long	80308	80092	80236	80020	80272	80056
	PLUS	Short	80324	80108	80252	80036	80288	80072
		Normal	80325	80109	80253	80037	80289	80073
		Long	80326	80110	80254	80038	80290	80074
XL	STD	Short	80309	80093	80237	80021	80273	80057
		Normal	80310	80094	80238	80022	80274	80058
		Long	80311	80095	80239	80023	80275	80059
	PLUS	Short	80327	80111	80255	80039	80291	80075
		Normal	80328	80112	80256	80040	80292	80076
		Long	80329	80113	80257	80041	80293	80077

●●● Closed toe ○○○ Open toe

Product codes

Tights (panty)




FOR THIS PRODUCT - Please add '-CM' to the end of your product codes when ordering

			Class 1						Class 2					
			Black		Caramel		Savannah		Black		Caramel		Savannah	
			●	○	●	○	●	○	●	○	●	○	●	○
XXS	STD	Normal	69297	69289	69265	69257	69329	69305	70857	70849	70825	70817	70889	70865
		Long	69298	69290	69266	69258	69330	69306	70858	70850	70826	70818	70890	70866
	PLUS	Normal	69301	69293	69269	69261	69341	69317	70861	70853	70829	70821	70901	70877
		Long	69302	69294	69270	69262	69342	69318	70862	70854	70830	70822	70902	70878
XS	STD	Normal	49433	49337	49401	49305	69331	69307	41064	41048	41000	40984	70891	70867
		Long	49434	49338	49402	49306	69332	69308	41065	41049	41001	40985	70892	70868
	PLUS	Normal	49441	49345	49409	49313	69343	69319	41072	41056	41008	40992	70903	70879
		Long	49442	49346	49410	49314	69344	69320	41073	41057	41009	40993	70904	70880
S	STD	Normal	49435	49339	49403	49307	69333	69309	41066	41050	41002	40986	70893	70869
		Long	49436	49340	49404	49308	69334	69310	41067	41051	41003	40987	70894	70870
	PLUS	Normal	49443	49347	49411	49315	69345	69321	41074	41058	41010	40994	70905	70881
		Long	49444	49348	49412	49316	69346	69322	41075	41059	41011	40995	70906	70882
M	STD	Normal	49437	49341	49405	49309	69335	69311	41068	41052	41004	40988	70895	70871
		Long	49438	49342	49406	49310	69336	69312	41069	41053	41005	40989	70896	70872
	PLUS	Normal	49445	49349	49413	49317	69347	69323	41076	41060	41012	40996	70907	70883
		Long	49446	49350	49414	49318	69348	69324	41077	41061	41013	40997	70908	70884
L	STD	Normal	49439	49343	49407	49311	69337	69313	41070	41054	41006	40990	70897	70873
		Long	49440	49344	49408	49312	69338	69314	41071	41055	41007	40991	70898	70874
	PLUS	Normal	49447	49351	49415	49319	69349	69325	41078	41062	41014	40998	70909	70885
		Long	49448	49352	49416	49320	69350	69326	41079	41063	41015	40999	70910	70886
XL	STD	Normal	69299	69291	69267	69259	69339	69315	70859	70851	70827	70819	70899	70875
		Long	69300	69292	69268	69260	69340	69316	70860	70852	70828	70820	70900	70876
	PLUS	Normal	69303	69295	69271	69263	69351	69327	70863	70855	70831	70823	70911	70887
		Long	69304	69296	69272	69264	69352	69328	70864	70856	70832	70824	70912	70888

●●● Closed toe ○○○ Open toe

Product codes

Tights (panty) extra




FOR THIS PRODUCT - Please add '-CM' to the end of your product codes when ordering

			Class 1						Class 2					
			Black		Caramel		Savannah		Black		Caramel		Savannah	
			●	○	●	○	●	○	●	○	●	○	●	○
XXS	STD	Normal	69409	69401	69377	69369	69441	69417	71049	71041	71017	71009	71081	71057
		Long	69410	69402	69378	69370	69442	69418	71050	71042	71018	71010	71082	71058
	PLUS	Normal	69413	69405	69381	69373	69453	69429	71053	71045	71021	71013	71093	71069
		Long	69414	69406	69382	69374	69454	69430	71054	71046	71022	71014	71094	71070
XS	STD	Normal	49817	49721	49785	49689	69443	69419	41448	41432	41384	41368	71083	71059
		Long	49818	49722	49786	49690	69444	69420	41449	41433	41385	41369	71084	71060
	PLUS	Normal	49825	49729	49793	49697	69455	69431	41456	41440	41392	41376	71095	71071
		Long	49826	49730	49794	49698	69456	69432	41457	41441	41393	41377	71096	71072
S	STD	Normal	49819	49723	49787	49691	69445	69421	41450	41434	41386	41370	71085	71061
		Long	49820	49724	49788	49692	69446	69422	41451	41435	41387	41371	71086	71062
	PLUS	Normal	49827	49731	49795	49699	69457	69433	41458	41442	41394	41378	71097	71073
		Long	49828	49732	49796	49700	69458	69434	41459	41443	41395	41379	71098	71074
M	STD	Normal	49821	49725	49789	49693	69447	69423	41452	41436	41388	41372	71087	71063
		Long	49822	49726	49790	49694	69448	69424	41453	41437	41389	41373	71088	71064
	PLUS	Normal	49829	49733	49797	49701	69459	69435	41460	41444	41396	41380	71099	71075
		Long	49830	49734	49798	49702	69460	69436	41461	41445	41397	41381	71100	71076
L	STD	Normal	49823	49727	49791	49695	69449	69425	41454	41438	41390	41374	71089	71065
		Long	49824	49728	49792	49696	69450	69426	41455	41439	41391	41375	71090	71066
	PLUS	Normal	49831	49735	49799	49703	69461	69437	41462	41446	41398	41382	71101	71077
		Long	49832	49736	49800	49704	69462	69438	41463	41447	41399	41383	71102	71078
XL	STD	Normal	69411	69403	69379	69371	69451	69427	71051	71043	71019	71011	71091	71067
		Long	69412	69404	69380	69372	69452	69428	71052	71044	71020	71012	71092	71068
	PLUS	Normal	69415	69407	69383	69375	69463	69439	71055	71047	71023	71015	71103	71079
		Long	69416	69408	69384	69376	69464	69440	71056	71048	71024	71016	71104	71080

●● Closed toe ○○ Open toe

Product codes

Tights (panty) - Bodyform



FOR THIS PRODUCT - Please add '-CM' to the end of your product codes when ordering

			Class 2					
			Black		Caramel		Savannah	
			●	○	●	○	●	○
XXS	STD	Normal	70937	70929	70921	70913	70969	70945
		Long	70938	70930	70922	70914	70970	70946
	PLUS	Normal	70941	70933	70925	70917	70981	70957
		Long	70942	70934	70926	70918	70982	70958
XS	STD	Normal	41256	41240	41192	41176	70971	70947
		Long	41257	41241	41193	41177	70972	70948
	PLUS	Normal	41264	41248	41200	41184	70983	70959
		Long	41265	41249	41201	41185	70984	70960
S	STD	Normal	41258	41242	41194	41178	70973	70949
		Long	41259	41243	41195	41179	70974	70950
	PLUS	Normal	41266	41250	41202	41186	70985	70961
		Long	41267	41251	41203	41187	70986	70962
M	STD	Normal	41260	41244	41196	41180	70975	70951
		Long	41261	41245	41197	41181	70976	70952
	PLUS	Normal	41268	41252	41204	41188	70987	70963
		Long	41269	41253	41205	41189	70988	70964
L	STD	Normal	41262	41246	41198	41182	70977	70953
		Long	41263	41247	41199	41183	70978	70954
	PLUS	Normal	41270	41254	41206	41190	70989	70965
		Long	41271	41255	41207	41191	70990	70966
XL	STD	Normal	70939	70931	70923	70915	70979	70955
		Long	70940	70932	70924	70916	70980	70956
	PLUS	Normal	70943	70935	70927	70919	70991	70967
		Long	70944	70936	70928	70920	70992	70968

●● Closed toe ○○ Open toe


Available in MTM

Offers tailored sizing, colour options, and style choices, ensuring the right product with the right fit.

Find out more in the Made to Measure section



Sensinnov grip top band on thigh



Knobbled grip top band

Material composition

All styles

52% Polyamide, 34% Elastane, 14% cotton (content can vary slightly by compression choice)

Latex free



Scan here for this product's measuring form

Essential THERMOREGULATING

With 100% cotton against the skin and properties that maintain an optimal skin climate for all year-round wear, this garment makes it ideal for wearers with venous or lymphatic conditions and sensitive skin.

Compression class

18-21 mmHg

RAL 1

23-32 mmHg

RAL 2

34-46 mmHg

RAL 3



Sensitive skin



100% cotton against the skin



Breathable



Available in MTM

Suitable for

- Venous Insufficiency; including leg ulceration, mixed disease
- Lymphoedema / chronic oedema
- Sensitive skin
- Lipoedema (lipalgia)
- Prevention of leg ulcers and oedema



Thigh with Sensinnov grip top in Black

Below knee in Natural with Open toe

Thigh with waist attachment in Natural

Tights (panty) in Black with Open toe

EMIS search tips

Product name

Essential THERMOREGULATING **SIG ESS THERMO**
EMIS TIP

Compression class (mmHg level)

Class 1 **18** Class 2 **23** Class 3 **34**
EMIS TIP EMIS TIP EMIS TIP

Sizing chart

All measurements are shown in cm


Below knee, thigh and tights (panty)

	XXS		XS		S		M		L		XL	
Ankle (cB)	18-20		20-22		22-24		24-26		26-29		29-32	
	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS
Calf (cC)	25-32	29-36	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49	42-49	43-51
Thigh (cG) with grip top	42-52	47-56	47-56	52-61	52-61	56-65	56-65	61-71	61-71	66-76	66-76	73-85
Thigh (cG) (For thigh w/ attachment and tights only)	42-54	47-59	47-59	52-64	52-64	56-68	56-68	61-74	61-74	66-79	66-79	76-85
Hip (cH) (For tights only)	75-90	80-95	80-95	85-100	85-100	95-110	95-110	105-120	105-120	115-130	115-130	125-140
Waist (cT) (For tights only)	60-75	65-80	65-80	75-90	75-90	85-100	85-100	95-110	95-110	105-120	105-120	115-130
EMIS TIP	SMA	SMA PLU	SMA	SMA PLU	SMA	SMA PLU	MED	MED PLU	LAR	LAR PLU	LAR	LAR PLU


Length

	Short	Normal	Long
Below knee (lD)	31-34	35-40	41-45
Thigh (lG)	56-64	65-74	75-84
Tights (panty) (lG)	56-64	65-74	75-84
EMIS TIP	SHO	NOR	LON


Below knee




Thigh



Thigh w/ attachment



Tights (panty)




Styles and colours

- Closed toe **CLO**
- Open toe **OPE**

		Class 1		Class 2		Class 3	
		Black	Nature	Black	Nature	Black	Nature
Below knee	BEL	●*○*	●*○*	●●	●●	●*○*	●*○*
Below knee with knobbled grip top	BEL GRI	-	-	●*○*	●*○*	-	-
Below knee maxi foot	BEL MAX	●*	-	●*	●*	●*	-
Thigh without grip top	THI	●●	●●	●●	●●	-	-
Thigh with knobbled grip top	THI GRI	●●	●●	●●	●●	●●	●●
Thigh with Sensinnov grip top	THI SEN	●●	●●	●●	●●	-	-
Thigh with waist attachment	THI WAI	●●	●●	●●	●●	●●	●●
Tights (panty)	PAN	●*○*	●*○*	●*○*	●●	●*○*	●●
EMIS TIP		BLK	NAT	BLK	NAT	BLK	NAT

* Not available in short length



Scan here for this product's ordering codes

Essential THERMOREGULATING



Product codes

Below knee

			Class 1				Class 2				Class 3			
			Black		Nature		Black		Nature		Black		Nature	
			●	○	●	○	●	○	●	○	●	○	●	○
XXS	STD	Short	-	-	-	-	63430	63406	63418	63394	-	-	-	-
		Normal	63971	63899	63935	56236	56298	56282	56290	56274	73931	73899	73923	73875
		Long	63972	63900	63936	56237	56299	56283	56291	56275	73932	73900	73924	73876
	PLUS	Short	-	-	-	-	63436	63412	63424	63400	-	-	-	-
		Normal	63989	63917	63953	63891	56302	56286	56294	56278	73935	73911	73927	73887
		Long	63990	63918	63954	63892	56303	56287	56295	56279	73936	73912	73928	73888
XS	STD	Short	-	-	-	-	63431	63407	63419	63395	-	-	-	-
		Normal	63974	63902	63938	54034	50366	50318	50350	50302	51054	73901	51038	73877
		Long	63975	63903	63939	54035	50367	50319	50351	50303	51055	73902	51039	73878
	PLUS	Short	-	-	-	-	63437	63413	63425	63401	-	-	-	-
		Normal	63992	63920	63956	54042	50374	50326	50358	50310	51062	73913	51046	73889
		Long	63993	63921	63957	54043	50375	50327	50359	50311	51063	73914	51047	73890
S	STD	Short	-	-	-	-	63432	63408	63420	63396	-	-	-	-
		Normal	63977	63905	63941	54036	50368	50320	50352	50304	51056	73903	51040	73879
		Long	63978	63906	63942	54037	50369	50321	50353	50305	51057	73904	51041	73880
	PLUS	Short	-	-	-	-	63438	63414	63426	63402	-	-	-	-
		Normal	63995	63923	63959	54044	50376	50328	50360	50312	51064	73915	51048	73891
		Long	63996	63924	63960	54045	50377	50329	50361	50313	51065	73916	51049	73892
M	STD	Short	-	-	-	-	63433	63409	63421	63397	-	-	-	-
		Normal	63980	63908	63944	54038	50370	50322	50354	50306	51058	73905	51042	73881
		Long	63981	63909	63945	54039	50371	50323	50355	50307	51059	73906	51043	73882
	PLUS	Short	-	-	-	-	63439	63415	63427	63403	-	-	-	-
		Normal	63998	63926	63962	54046	50378	50330	50362	50314	51066	73917	51050	73893
		Long	63999	63927	63963	54047	50379	50331	50363	50315	51067	73918	51051	73894
L	STD	Short	-	-	-	-	63434	63410	63422	63398	-	-	-	-
		Normal	63983	63911	63947	54040	50372	50324	50356	50308	51060	73907	51044	73883
		Long	63984	63912	63948	54041	50373	50325	50357	50309	51061	73908	51045	73884
	PLUS	Short	-	-	-	-	63440	63416	63428	63404	-	-	-	-
		Normal	64001	63929	63965	54048	50380	50332	50364	50316	51068	73919	51052	73895
		Long	64002	63930	63966	54049	50381	50333	50365	50317	51069	73920	51053	73896
XL	STD	Short	-	-	-	-	63435	63411	63423	63399	-	-	-	-
		Normal	63986	63914	63950	56238	56300	56284	56292	56276	73933	73909	73925	73885
		Long	63987	63915	63951	56239	56301	56285	56293	56277	73934	73910	73926	73886
	PLUS	Short	-	-	-	-	63441	63417	63429	63405	-	-	-	-
		Normal	64004	63932	63968	56240	56304	56288	56296	56280	73937	73921	73929	73897
		Long	64005	63933	63969	56241	56305	56289	56297	56281	73938	73922	73930	73898

● Closed toe ○ Open toe

Any questions or queries? Contact our Customer Service team on 01264 326 666 or email customerservices@sigvaris.com

Product codes

Below knee with knobbed grip top

			Class 2				
			Black		Nature		
			●	○	●	○	
XXS	STD	Normal	64475	64453	56264	56262	
		Long	64476	64454	56265	56263	
		PLUS	Normal	64666	64662	64664	64660
	XS	STD	Normal	64477	64455	50398	50382
			Long	64478	64456	50399	50383
			PLUS	Normal	64487	64465	50406
S	STD	Normal	64479	64457	50400	50384	
		Long	64480	64458	50401	50385	
		PLUS	Normal	64489	64467	50408	50392
M	STD	Normal	64481	64459	50402	50386	
		Long	64482	64460	50403	50387	
		PLUS	Normal	64491	64469	50410	50394
L	STD	Normal	64483	64461	50404	50388	
		Long	64484	64462	50405	50389	
		PLUS	Normal	64493	64471	50412	50396
XL	STD	Normal	64485	64463	56270	56266	
		Long	64486	64464	56271	56267	
		PLUS	Normal	64495	64473	56272	56268
XL	STD	Long	64496	64474	56273	56269	

● Closed toe ○ Open toe

Below knee maxi foot (Shoe size 10-13)


			Class 1	Class 2		Class 3		
			Black	Black	Nature	Black		
			●	●	●	●		
XS	STD	Normal	50862	50990	50974	51118		
		Long	50863	50991	50975	51119		
		PLUS	Normal	50870	50998	50982	51126	
	S	STD	Long	50871	50999	50983	51127	
			PLUS	Normal	50864	50992	50976	51120
				Long	50865	50993	50977	51121
M	STD	Normal		50872	51000	50984	51128	
		Long	50873	51001	50985	51129		
		L	STD	Normal	50866	50994	50978	51122
Long	50867			50995	50979	51123		
XL	STD			Normal	50874	51002	50986	51130
		Long	50875	51003	50987	51131		
		XL	STD	Normal	50868	50996	50980	51124
Long	50869			50997	50981	51125		
PLUS	Normal			50876	51004	50988	51132	
	Long		50877	51005	50989	51133		
	PLUS		Normal	56258	56322	56318	-	
Long			56259	56323	56319	-		
Long		56260	56324	56320	-			
XL	STD	Long	56261	56325	56321	-		

● Closed toe

Any questions or queries? Contact our Customer Service team on 01264 326 666 or email customerservices@sigvaris.com

Product codes

Thigh without grip top




			Class 1				Class 2			
			Black		Nature		Black		Nature	
			●	○	●	○	●	○	●	○
XXS	STD	Short	64114	64042	64078	64006	64605	64533	64569	64497
		Normal	64115	64043	64079	64007	64606	64534	64570	64498
		Long	64116	64044	64080	64008	64607	64535	64571	64499
	PLUS	Short	64132	64060	64096	64024	64623	64551	64587	64515
		Normal	64133	64061	64097	64025	64624	64552	64588	64516
		Long	64134	64062	64098	64026	64625	64553	64589	64517
XS	STD	Short	64117	64045	64081	64009	64608	64536	64572	64500
		Normal	64118	64046	64082	64010	64609	64537	64573	64501
		Long	64119	64047	64083	64011	64610	64538	64574	64502
	PLUS	Short	64135	64063	64099	64027	64626	64554	64590	64518
		Normal	64136	64064	64100	64028	64627	64555	64591	64519
		Long	64137	64065	64101	64029	64628	64556	64592	64520
S	STD	Short	64120	64048	64084	64012	64611	64539	64575	64503
		Normal	64121	64049	64085	64013	64612	64540	64576	64504
		Long	64122	64050	64086	64014	64613	64541	64577	64505
	PLUS	Short	64138	64066	64102	64030	64629	64557	64593	64521
		Normal	64139	64067	64103	64031	64630	64558	64594	64522
		Long	64140	64068	64104	64032	64631	64559	64595	64523
M	STD	Short	64123	64051	64087	64015	64614	64542	64578	64506
		Normal	64124	64052	64088	64016	64615	64543	64579	64507
		Long	64125	64053	64089	64017	64616	64544	64580	64508
	PLUS	Short	64141	64069	64105	64033	64632	64560	64596	64524
		Normal	64142	64070	64106	64034	64633	64561	64597	64525
		Long	64143	64071	64107	64035	64634	64562	64598	64526
L	STD	Short	64126	64054	64090	64018	64617	64545	64581	64509
		Normal	64127	64055	64091	64019	64618	64546	64582	64510
		Long	64128	64056	64092	64020	64619	64547	64583	64511
	PLUS	Short	64144	64072	64108	64036	64635	64563	64599	64527
		Normal	64145	64073	64109	64037	64636	64564	64600	64528
		Long	64146	64074	64110	64038	64637	64565	64601	64529
XL	STD	Short	64129	64057	64093	64021	64620	64548	64584	64512
		Normal	64130	64058	64094	64022	64621	64549	64585	64513
		Long	64131	64059	64095	64023	64622	64550	64586	64514
	PLUS	Short	64147	64075	64111	64039	64638	64566	64602	64530
		Normal	64148	64076	64112	64040	64639	64567	64603	64531
		Long	64149	64077	64113	64041	64640	64568	64604	64532

● Closed toe ○ Open toe

Product codes

Thigh with knobbed grip top




			Class 1				Class 2				Class 3				
			Black		Nature		Black		Nature		Black		Nature		
			●	○	●	○	●	○	●	○	●	○	●	○	
XXS	STD	Short	64225	64153	64189	56326	56371	56347	56359	56335	76661	76625	74052	74016	
		Normal	64226	64154	64190	56327	56372	56348	56360	56336	76662	76626	74053	74017	
		Long	64227	64155	64191	56328	56373	56349	56361	56337	76663	76627	74054	74018	
	PLUS	Short	64243	64171	64207	64150	56377	56353	56365	56341	76679	76643	74070	74034	
		Normal	64244	64172	64208	64151	56378	56354	56366	56342	76680	76644	74071	74035	
		Long	64245	64173	64209	64152	56379	56355	56367	56343	76681	76645	74072	74036	
	XS	STD	Short	64228	64156	64192	50262	50534	50462	50510	50438	76664	76628	74055	74019
			Normal	64229	64157	64193	50263	50535	50463	50511	50439	76665	76629	74056	74020
			Long	64230	64158	64194	50264	50536	50464	50512	50440	76666	76630	74057	74021
PLUS		Short	64246	64174	64210	50274	50546	50474	50522	50450	76682	76646	74073	74037	
		Normal	64247	64175	64211	50275	50547	50475	50523	50451	76683	76647	74074	74038	
		Long	64248	64176	64212	50276	50548	50476	50524	50452	76684	76648	74075	74039	
S	STD	Short	64231	64159	64195	50265	50537	50465	50513	50441	76667	76631	74058	74022	
		Normal	64232	64160	64196	50266	50538	50466	50514	50442	76668	76632	74059	74023	
		Long	64233	64161	64197	50267	50539	50467	50515	50443	76669	76633	74060	74024	
	PLUS	Short	64249	64177	64213	50277	50549	50477	50525	50453	76685	76649	74076	74040	
		Normal	64250	64178	64214	50278	50550	50478	50526	50454	76686	76650	74077	74041	
		Long	64251	64179	64215	50279	50551	50479	50527	50455	76687	76651	74078	74042	
M	STD	Short	64234	64162	64198	50268	50540	50468	50516	50444	76670	76634	74061	74025	
		Normal	64235	64163	64199	50269	50541	50469	50517	50445	76671	76635	74062	74026	
		Long	64236	64164	64200	50270	50542	50470	50518	50446	76672	76636	74063	74027	
	PLUS	Short	64252	64180	64216	50280	50552	50480	50528	50456	76688	76652	74079	74043	
		Normal	64253	64181	64217	50281	50553	50481	50529	50457	76689	76653	74080	74044	
		Long	64254	64182	64218	50282	50554	50482	50530	50458	76690	76654	74081	74045	
L	STD	Short	64237	64165	64201	50271	50543	50471	50519	50447	76673	76637	74064	74028	
		Normal	64238	64166	64202	50272	50544	50472	50520	50448	76674	76638	74065	74029	
		Long	64239	64167	64203	50273	50545	50473	50521	50449	76675	76639	74066	74030	
	PLUS	Short	64255	64183	64219	50283	50555	50483	50531	50459	76691	76655	74082	74046	
		Normal	64256	64184	64220	50284	50556	50484	50532	50460	76692	76656	74083	74047	
		Long	64257	64185	64221	50285	50557	50485	50533	50461	76693	76657	74084	74048	
XL	STD	Short	64240	64168	64204	56329	56374	56350	56362	56338	76676	76640	74067	74031	
		Normal	64241	64169	64205	56330	56375	56351	56363	56339	76677	76641	74068	74032	
		Long	64242	64170	64206	56331	56376	56352	56364	56340	76678	76642	74069	74033	
	PLUS	Short	64258	64186	64222	56332	56380	56356	56368	56344	76694	76658	74085	74049	
		Normal	64259	64187	64223	56333	56381	56357	56369	56345	76695	76659	74086	74050	
		Long	64260	64188	64224	56334	56382	56358	56370	56346	76696	76660	74087	74051	

● Closed toe ○ Open toe

Product codes

Thigh with seninnov grip top




			Class 1				Class 2			
			Black		Nature		Black		Nature	
			●	○	●	○	●	○	●	○
XXS	STD	Short	73119	73047	73083	73011	73263	73191	73227	73155
		Normal	73120	73048	73084	73012	73264	73192	73228	73156
		Long	73121	73049	73085	73013	73265	73193	73229	73157
	PLUS	Short	73137	73065	73101	73029	73281	73209	73245	73173
		Normal	73138	73066	73102	73030	73282	73210	73246	73174
		Long	73139	73067	73103	73031	73283	73211	73247	73175
XS	STD	Short	73122	73050	73086	73014	73266	73194	73230	73158
		Normal	73123	73051	73087	73015	73267	73195	73231	73159
		Long	73124	73052	73088	73016	73268	73196	73232	73160
	PLUS	Short	73140	73068	73104	73032	73284	73212	73248	73176
		Normal	73141	73069	73105	73033	73285	73213	73249	73177
		Long	73142	73070	73106	73034	73286	73214	73250	73178
S	STD	Short	73125	73053	73089	73017	73269	73197	73233	73161
		Normal	73126	73054	73090	73018	73270	73198	73234	73162
		Long	73127	73055	73091	73019	73271	73199	73235	73163
	PLUS	Short	73143	73071	73107	73035	73287	73215	73251	73179
		Normal	73144	73072	73108	73036	73288	73216	73252	73180
		Long	73145	73073	73109	73037	73289	73217	73253	73181
M	STD	Short	73128	73056	73092	73020	73272	73200	73236	73164
		Normal	73129	73057	73093	73021	73273	73201	73237	73165
		Long	73130	73058	73094	73022	73274	73202	73238	73166
	PLUS	Short	73146	73074	73110	73038	73290	73218	73254	73182
		Normal	73147	73075	73111	73039	73291	73219	73255	73183
		Long	73148	73076	73112	73040	73292	73220	73256	73184
L	STD	Short	73131	73059	73095	73023	73275	73203	73239	73167
		Normal	73132	73060	73096	73024	73276	73204	73240	73168
		Long	73133	73061	73097	73025	73277	73205	73241	73169
	PLUS	Short	73149	73077	73113	73041	73293	73221	73257	73185
		Normal	73150	73078	73114	73042	73294	73222	73258	73186
		Long	73151	73079	73115	73043	73295	73223	73259	73187
XL	STD	Short	73134	73062	73098	73026	73278	73206	73242	73170
		Normal	73135	73063	73099	73027	73279	73207	73243	73171
		Long	73136	73064	73100	73028	73280	73208	73244	73172
	PLUS	Short	73152	73080	73116	73044	73296	73224	73260	73188
		Normal	73153	73081	73117	73045	73297	73225	73261	73189
		Long	73154	73082	73118	73046	73298	73226	73262	73190

● Closed toe ○ Open toe

Product codes

Thigh with waist attachment - Left




			Class 1				Class 2				Class 3				
			Black		Nature		Black		Nature		Black		Nature		
			●	○	●	○	●	○	●	○	●	○	●	○	
XXS	STD	Short	78494	78350	78458	78314	78782	78638	78746	78602	79070	78926	79034	78890	
		Normal	78495	78351	78459	78315	78783	78639	78747	78603	79071	78927	79035	78891	
		Long	78496	78352	78460	78316	78784	78640	78748	78604	79072	78928	79036	78892	
	PLUS	Short	78512	78368	78476	78332	78800	78656	78764	78620	79088	78944	79052	78908	
		Normal	78513	78369	78477	78333	78801	78657	78765	78621	79089	78945	79053	78909	
		Long	78514	78370	78478	78334	78802	78658	78766	78622	79090	78946	79054	78910	
	XS	STD	Short	78497	78353	78461	78317	78785	78641	78749	78605	79073	78929	79037	78893
			Normal	78498	78354	78462	78318	78786	78642	78750	78606	79074	78930	79038	78894
			Long	78499	78355	78463	78319	78787	78643	78751	78607	79075	78931	79039	78895
PLUS		Short	78515	78371	78479	78335	78803	78659	78767	78623	79091	78947	79055	78911	
		Normal	78516	78372	78480	78336	78804	78660	78768	78624	79092	78948	79056	78912	
		Long	78517	78373	78481	78337	78805	78661	78769	78625	79093	78949	79057	78913	
S	STD	Short	78500	78356	78464	78320	78788	78644	78752	78608	79076	78932	79040	78896	
		Normal	78501	78357	78465	78321	78789	78645	78753	78609	79077	78933	79041	78897	
		Long	78502	78358	78466	78322	78790	78646	78754	78610	79078	78934	79042	78898	
	PLUS	Short	78518	78374	78482	78338	78806	78662	78770	78626	79094	78950	79058	78914	
		Normal	78519	78375	78483	78339	78807	78663	78771	78627	79095	78951	79059	78915	
		Long	78520	78376	78484	78340	78808	78664	78772	78628	79096	78952	79060	78916	
M	STD	Short	78503	78359	78467	78323	78791	78647	78755	78611	79079	78935	79043	78899	
		Normal	78504	78360	78468	78324	78792	78648	78756	78612	79080	78936	79044	78900	
		Long	78505	78361	78469	78325	78793	78649	78757	78613	79081	78937	79045	78901	
	PLUS	Short	78521	78377	78485	78341	78809	78665	78773	78629	79097	78953	79061	78917	
		Normal	78522	78378	78486	78342	78810	78666	78774	78630	79098	78954	79062	78918	
		Long	78523	78379	78487	78343	78811	78667	78775	78631	79099	78955	79063	78919	
L	STD	Short	78506	78362	78470	78326	78794	78650	78758	78614	79082	78938	79046	78902	
		Normal	78507	78363	78471	78327	78795	78651	78759	78615	79083	78939	79047	78903	
		Long	78508	78364	78472	78328	78796	78652	78760	78616	79084	78940	79048	78904	
	PLUS	Short	78524	78380	78488	78344	78812	78668	78776	78632	79100	78956	79064	78920	
		Normal	78525	78381	78489	78345	78813	78669	78777	78633	79101	78957	79065	78921	
		Long	78526	78382	78490	78346	78814	78670	78778	78634	79102	78958	79066	78922	
XL	STD	Short	78509	78365	78473	78329	78797	78653	78761	78617	79085	78941	79049	78905	
		Normal	78510	78366	78474	78330	78798	78654	78762	78618	79086	78942	79050	78906	
		Long	78511	78367	78475	78331	78799	78655	78763	78619	79087	78943	79051	78907	
	PLUS	Short	78527	78383	78491	78347	78815	78671	78779	78635	79103	78959	79067	78923	
		Normal	78528	78384	78492	78348	78816	78672	78780	78636	79104	78960	79068	78924	
		Long	78529	78385	78493	78349	78817	78673	78781	78637	79105	78961	79069	78925	

● Closed toe ○ Open toe

Product codes

Thigh with waist attachment - Right




			Class 1				Class 2				Class 3			
			Black		Nature		Black		Nature		Black		Nature	
			●	○	●	○	●	○	●	○	●	○	●	○
XXS	STD	Short	78566	78422	78530	78386	78854	78710	78818	78674	79142	78998	79106	78962
		Normal	78567	78423	78531	78387	78855	78711	78819	78675	79143	78999	79107	78963
		Long	78568	78424	78532	78388	78856	78712	78820	78676	79144	79000	79108	78964
	PLUS	Short	78584	78440	78548	78404	78872	78728	78836	78692	79160	79016	79124	78980
		Normal	78585	78441	78549	78405	78873	78729	78837	78693	79161	79017	79125	78981
		Long	78586	78442	78550	78406	78874	78730	78838	78694	79162	79018	79126	78982
XS	STD	Short	78569	78425	78533	78389	78857	78713	78821	78677	79145	79001	79109	78965
		Normal	78570	78426	78534	78390	78858	78714	78822	78678	79146	79002	79110	78966
		Long	78571	78427	78535	78391	78859	78715	78823	78679	79147	79003	79111	78967
	PLUS	Short	78587	78443	78551	78407	78875	78731	78839	78695	79163	79019	79127	78983
		Normal	78588	78444	78552	78408	78876	78732	78840	78696	79164	79020	79128	78984
		Long	78589	78445	78553	78409	78877	78733	78841	78697	79165	79021	79129	78985
S	STD	Short	78572	78428	78536	78392	78860	78716	78824	78680	79148	79004	79112	78968
		Normal	78573	78429	78537	78393	78861	78717	78825	78681	79149	79005	79113	78969
		Long	78574	78430	78538	78394	78862	78718	78826	78682	79150	79006	79114	78970
	PLUS	Short	78590	78446	78554	78410	78878	78734	78842	78698	79166	79022	79130	78986
		Normal	78591	78447	78555	78411	78879	78735	78843	78699	79167	79023	79131	78987
		Long	78592	78448	78556	78412	78880	78736	78844	78700	79168	79024	79132	78988
M	STD	Short	78575	78431	78539	78395	78863	78719	78827	78683	79151	79007	79115	78971
		Normal	78576	78432	78540	78396	78864	78720	78828	78684	79152	79008	79116	78972
		Long	78577	78433	78541	78397	78865	78721	78829	78685	79153	79009	79117	78973
	PLUS	Short	78593	78449	78557	78413	78881	78737	78845	78701	79169	79025	79133	78989
		Normal	78594	78450	78558	78414	78882	78738	78846	78702	79170	79026	79134	78990
		Long	78595	78451	78559	78415	78883	78739	78847	78703	79171	79027	79135	78991
L	STD	Short	78578	78434	78542	78398	78866	78722	78830	78686	79154	79010	79118	78974
		Normal	78579	78435	78543	78399	78867	78723	78831	78687	79155	79011	79119	78975
		Long	78580	78436	78544	78400	78868	78724	78832	78688	79156	79012	79120	78976
	PLUS	Short	78596	78452	78560	78416	78884	78740	78848	78704	79172	79028	79136	78992
		Normal	78597	78453	78561	78417	78885	78741	78849	78705	79173	79029	79137	78993
		Long	78598	78454	78562	78418	78886	78742	78850	78706	79174	79030	79138	78994
XL	STD	Short	78581	78437	78545	78401	78869	78725	78833	78689	79157	79013	79121	78977
		Normal	78582	78438	78546	78402	78870	78726	78834	78690	79158	79014	79122	78978
		Long	78583	78439	78547	78403	78871	78727	78835	78691	79159	79015	79123	78979
	PLUS	Short	78599	78455	78563	78419	78887	78743	78851	78707	79175	79031	79139	78995
		Normal	78600	78456	78564	78420	78888	78744	78852	78708	79176	79032	79140	78996
		Long	78601	78457	78565	78421	78889	78745	78853	78709	79177	79033	79141	78997

● Closed toe ○ Open toe

Product codes

Tights (panty)



			Class 1				Class 2				Class 3			
			Black		Nature		Black		Nature		Black		Nature	
			●	○	●	○	●	○	●	○	●	○	●	○
XXS	STD	Short	-	-	-	-	-	-	74004	73992	-	-	74124	74088
		Normal	64333	64285	64309	64261	56407	56391	56399	56383	76721	76697	74125	74089
		Long	64334	64286	64310	64262	56408	56392	56400	56384	76722	76698	74126	74090
	PLUS	Short	-	-	-	-	-	-	74010	73998	-	-	74142	74106
		Normal	64345	64297	64321	64273	56411	56395	56403	56387	76733	76709	74143	74107
		Long	64346	64298	64322	64274	56412	56396	56404	56388	76734	76710	74144	74108
XS	STD	Short	-	-	-	-	-	-	74005	73993	-	-	74127	74091
		Normal	64335	64287	64311	64263	50638	50590	50622	50574	76723	76699	74128	74092
		Long	64336	64288	64312	64264	50639	50591	50623	50575	76724	76700	74129	74093
	PLUS	Short	-	-	-	-	-	-	74011	73999	-	-	74145	74109
		Normal	64347	64299	64323	64275	50646	50598	50630	50582	76735	76711	74146	74110
		Long	64348	64300	64324	64276	50647	50599	50631	50583	76736	76712	74147	74111
S	STD	Short	-	-	-	-	-	-	74006	73994	-	-	74130	74094
		Normal	64337	64289	64313	64265	50640	50592	50624	50576	76725	76701	74131	74095
		Long	64338	64290	64314	64266	50641	50593	50625	50577	76726	76702	74132	74096
	PLUS	Short	-	-	-	-	-	-	74012	74000	-	-	74148	74112
		Normal	64349	64301	64325	64277	50648	50600	50632	50584	76737	76713	74149	74113
		Long	64350	64302	64326	64278	50649	50601	50633	50585	76738	76714	74150	74114
M	STD	Short	-	-	-	-	-	-	74007	73995	-	-	74133	74097
		Normal	64339	64291	64315	64267	50642	50594	50626	50578	76727	76703	74134	74098
		Long	64340	64292	64316	64268	50643	50595	50627	50579	76728	76704	74135	74099
	PLUS	Short	-	-	-	-	-	-	74013	74001	-	-	74151	74115
		Normal	64351	64303	64327	64279	50650	50602	50634	50586	76739	76715	74152	74116
		Long	64352	64304	64328	64280	50651	50603	50635	50587	76740	76716	74153	74117
L	STD	Short	-	-	-	-	-	-	74008	73996	-	-	74136	74100
		Normal	64341	64293	64317	64269	50644	50596	50628	50580	76729	76705	74137	74101
		Long	64342	64294	64318	64270	50645	50597	50629	50581	76730	76706	74138	74102
	PLUS	Short	-	-	-	-	-	-	74014	74002	-	-	74154	74118
		Normal	64353	64305	64329	64281	50652	50604	50636	50588	76741	76717	74155	74119
		Long	64354	64306	64330	64282	50653	50605	50637	50589	76742	76718	74156	74120
XL	STD	Short	-	-	-	-	-	-	74009	73997	-	-	74139	74103
		Normal	64343	64295	64319	64271	56409	56393	56401	56385	76731	76707	74140	74104
		Long	64344	64296	64320	64272	56410	56394	56402	56386	76732	76708	74141	74105
	PLUS	Short	-	-	-	-	-	-	74015	74003	-	-	74157	74121
		Normal	64355	64307	64331	64283	56413	56397	56405	56389	76743	76719	74158	74122
		Long	64356	64308	64332	64284	56414	56398	56406	56390	76744	76720	74159	74123

● Closed toe ○ Open toe



Material composition

75% Polyamide, 14%
Aramid Kermel, 11%
Elastane
Latex free

*Black /
Turquoise*

Black / Red



Scan here
for this product's
measuring form

Active WORK WEAR

These durable socks have been specially designed to protect feet in safety shoes, with highly resistant fabric to support the toe and heel and reinforced sole to add extra protection to the ball of the foot.

Compression class

15-20
mmHg

AFNOR
2

Product benefits



Foot shell
protection Reinforced
zone Anti-
bacterial

Suitable for

- Venous insufficiency
- Mild lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Mixed disease

EMIS search tips

Product name
Active WORK WEAR

SIG ACT WORK
EMIS TIP

Sizing chart

All measurements are shown in cm

Below knee


	XS*	S	M	L	XL
Ankle (cB)	17-19	19-22	22-24	24-26	26-29
Calf (cC)	27-35	30-38	33-41	37-45	40-48
EMIS TIP	SMA	SMA	MED	LAR	LAR

*XS size only available in Blue / Turquoise

Length


		Normal	Long
Below knee (lD)		≤ 40	≥ 40
Shoe size	Female	36/40 - 3.5/6.5	37/41 - 4/7.5
	Male	40/43- 6/9	43/46- 9/11
EMIS TIP		NOR	LON

Below knee



Styles and colours

- Closed toe



		Class 2	
		Black / Turquoise	Black / Red
Below knee	FEM	●	-
	MAL	-	●



Scan here
for this
product's
ordering
codes

Active WORK WEAR




Product codes


Below knee



Female

		XS		S		M		L		XL	
		Normal	Long	Normal	Long	Normal	Long	Normal	Long	Normal	Long
	●	307239	307244	307240	307245	307241	307246	307242	307247	307243	307248

Male

		S		M		L		XL	
		Normal	Long	Normal	Long	Normal	Long	Normal	Long
	●	307249	307253	307250	307254	307251	307255	307252	307256

Light Grey

Material composition
 54% Polyamide, 32%
 Bamboo, 11% Elastane,
 3% Other fibres
 Latex free



Scan here
 for this product's
 measuring form

Active LEISURE

These socks combine a patented soft-knit bamboo fabric with ventilated calves and anti-bacterial properties to deliver long-lasting softness and comfort.

Compression class

15-20
mmHg

AFNOR
2

Product benefits



Suitable for

- Venous insufficiency
- Mild lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Mixed disease



Denim



Dusty Pink



Dark Grey

EMIS search tips

Product name
 Active WORK WEAR

SIG ACT WORK
 EMIS TIP

Sizing chart

All measurements are shown in cm

Female

Below knee

	XS	S	M	L	XL	XXL
Ankle (cB)	17-19	19-22	22-24	24-26	26-29	29-32
Calf (cC)	27-35	30-38	33-41	37-45	40-48	43-51
EMIS TIP	SMA	SMA	MED	LAR	LAR	LAR

Length

	Normal	Long
Below knee (lD)	≤ 40	≥ 40
Shoe size	36/40 - 3.5/6.5	37/41 - 4/7.5
EMIS TIP	NOR	LON

Male


Below knee

	S	M	L	XL	XXL
Ankle (cB)	19-22	22-24	24-26	26-29	29-32
Calf (cC)	31-38	34-41	36-43	39-46	42-49
EMIS TIP	SMA	MED	LAR	LAR	LAR

Length

	Normal	Long
Below knee (lD)	≤ 42	≥ 42
Shoe size	40/43 - 6/9	43/46 - 9/11
EMIS TIP	NOR	LON

Below knee



Styles and colours

●●●● Closed toe

		Class 2				
		Black	Denim	Dark Grey	Light Grey	Dusty Pink
Below knee	EMIS TIP					
	FEM	●	●	-	●	●
	MAL	●	-	●	-	-




Scan here
 for this
 product's
 ordering
 codes

Active LEISURE



Product codes

Below knee



Female

		XS		S		M		L		XL		XXL	
		Normal	Long	Normal	Long	Normal	Long	Normal	Long	Normal	Long	Normal	Long
Black	●	309171	309177	309172	309178	309173	309179	309174	309180	309175	309181	309176	309182
Light Grey	●	309183	309189	309184	309190	309185	309191	309186	309192	309187	309193	309188	309194
Denim	●	309195	309201	309196	309202	309197	309203	309198	309204	309199	309205	309200	309206
Light Pink	●	309207	309213	309208	309214	309209	309215	309210	309216	309211	309217	309212	309218

●●●● Closed toe

Male

		S		M		L		XL		XXL	
		Normal	Long	Normal	Long	Normal	Long	Normal	Long	Normal	Long
Black	●	309137	309142	309138	309143	309139	309144	309140	309145	309141	309146
Grey	●	309147	309156	309152	309157	309153	309158	309154	309159	309155	309160

Available in MTM

Offers tailored sizing, colour options, and style choices, ensuring the right product with the right fit.

Find out more in the Made to Measure section



Knobbed grip
top band on thigh

Thigh with knobbed grip top



Below knee



Thigh with waist attachment

Material composition

All styles

70% Polyamide, 30% Natural rubber latex



Scan here for this product's measuring form

TRADITIONAL

Providing high-containment for wearers with more severe conditions, its durable fabric and consistent pressure make it a reliable choice for intensive, long-term therapy in clinical practice.

Compression class

23-32 mmHg 34-46 mmHg

RAL 2 RAL 3

Product benefits



Natural rubber latex Durable Containment Made to Measure

Suitable for

- Lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Lipoedema (lipalgia)
- Active leg ulcers

EMIS search tips

Product name
TRADITIONAL

SIG TRA
EMIS TIP

Compression class (mmHg level)
Class 2 23 **Class 3 34**
EMIS TIP EMIS TIP

Sizing chart

All measurements are shown in cm

Below knee, thigh and thigh with waist attachment


	XS	S		M		L		XL	
Ankle (cB)	17-19	19-23		23-26		26-29		29-32	
	STD	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS
Calf (cC)	28-31	31-36	35-40	34-39	39-44	37-42	43-48	40-45	47-52
Thigh (cG) with grip top	40-50	42-52	50-60	50-60	58-68	58-68	66-76	-	-
Thigh (cG) (For thigh w/ attachment only)	40-50*	42-52	50-60	50-60	58-68	58-68	66-76	61-71*	69-79*
Hip (cH)	-	up to 110	up to 110	up to 120	up to 120	up to 130	up to 130	-	-
Waist (cT)	-	up to 80*	up to 80	up to 90*	up to 90	up to 100*	up to 100	-	-
EMIS TIP	SMA	SMA	SMA PLU	MED	MED PLU	LAR	LAR PLU	LAR	LAR PLU

* Sizes only available in Compression Class 2


Length

	Short	Long
Below knee (ℓD)	33-37	38-44
Thigh (ℓG)	61-70	71-80
Thigh with waist attachment (ℓG)	61-70	71-80
EMIS TIP	SHO	LON


Below knee




Thigh



Thigh w/ attachment



Tights (panty)



Styles and colours

	EMIS TIP	Class 2	Class 3
		Beige	Beige
Below knee	BEL	○	○
Thigh without grip top	THI	○	○
Thigh with knobbed grip top	THI GRI	○	○
Thigh with waist attachment (Left and right)	THI WAI LEF	○	○
	THI WAI RIG	○	○

○ Open toe



Scan here for this product's ordering codes

TRADITIONAL




Product codes

Below knee



			Class 2	Class 3
			Beige	Beige
			○	○
XS	STD	Short	10103	-
		Long	10104	-
S	STD	Short	10105	10271
		Long	10106	10272
	PLUS	Short	10114	10280
		Long	10115	10281
M	STD	Short	10107	10273
		Long	10108	10274
	PLUS	Short	10116	10282
		Long	10117	10283
L	STD	Short	10109	10275
		Long	10110	10276
	PLUS	Short	10118	10284
		Long	10119	10285
XL	STD	Short	10111	10277
		Long	10112	10278
	PLUS	Short	10120	10286
		Long	10121	10287

Thigh without grip top




			Class 2	Class 3
			Beige	Beige
			○	○
XS	STD	Long	10148	-
		Short	15168	15174
S	STD	Long	10149	10312
		Short	15171	15177
	PLUS	Short	15171	15177
		Long	10154	10317
M	STD	Short	15169	15175
		Long	10150	10313
	PLUS	Short	15172	15178
		Long	10155	10318
L	STD	Short	15170	15176
		Long	10151	10314
	PLUS	Short	15173	15179
		Long	10156	10319
XL	STD	Long	10152	-
	PLUS	Long	10157	-

○ Open toe

○ Open toe

Product codes


Thigh with grip top



			Class 2	Class 3
			Beige	Beige
			○	○
S	STD	Short	18833	18839
		Long	15112	15132
	PLUS	Short	18836	18842
		Long	15116	15136
M	STD	Short	18834	18840
		Long	15113	15133
	PLUS	Short	18837	18843
		Long	15117	15137
L	STD	Short	18835	18841
		Long	15114	15134
	PLUS	Short	18838	18844
		Long	15118	15138

○ Open toe

Thigh with waist attachment



			Left		Right	
			Class 2	Class 3	Class 2	Class 3
			Beige	Beige	Beige	Beige
			○	○	○	○
S	STD	Short	80810	-	80822	-
		Long	80811	-	80823	-
	PLUS	Short	80816	80834	80828	80840
		Long	80817	80835	80829	80841
M	STD	Short	80812	-	80824	-
		Long	80813	-	80825	-
	PLUS	Short	80818	80836	80830	80842
		Long	80819	80837	80831	80843
L	STD	Short	80814	-	80826	-
		Long	80815	-	80827	-
	PLUS	Short	80820	80838	80832	80844
		Long	80821	80839	80833	80845

○ Open toe



Material composition

Sleeve

70% Polyamide, 30% Natural rubber latex

Liner

62% Polyamide, 26% Elastane, 12% Cotton



Scan here for this product's measuring form

ULCER X Kit


The combination of compression sleeve and liner supports quicker healing and high compression levels in a comfortable and convenient manner.¹

Compression levels

$$23-32 \text{ mmHg} + 15-20 \text{ mmHg} = >45 \text{ mmHg}$$

Compression sleeve (1 device) Liner (1 pair) Kit

Product benefits



Suitable for

- Active venous ulcers
- Venous insufficiency

EMIS search tips

Product name
ULCER X Kit

SIG ULCER
EMIS TIP

Sizing chart

All measurements are shown in cm


Below knee

	S		M		L		XL	
Ankle (cB)	19-23		23-26		26-29		29-32	
Calf (cC)	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS
	31-36	35-40	34-39	39-44	37-42	43-48	40-45	47-52
EMIS TIP	SMA	SMA	MED	MED	LAR	LAR	LAR	LAR

Length

	Short	Long
Below knee (lD)	33-37	38-44
EMIS TIP	NOR	LON

Below knee



Styles and colours

- Closed toe



Beige

	EMIS TIP	
ULCER X Kit	STO	●



Scan here for this product's ordering codes



Material composition

Liner

62% Polyamide, 26% Elastane, 12% Cotton



Scan here for this product's measuring form

ULCER X Liner

The liner has a cotton blend that provides comfort against the skin and offers 15-20mmHg of compression.

Compression class

15-20 mmHg

Liner

Product benefits



Contains cotton Soft Breathable

Suitable for

- Early intervention of leg ulcers
- Venous conditions
- Prevention of ulcer recurrence
- Active venous ulcers

EMIS search tips

Product name

ULCER X Liner

SIG ULCER

EMIS TIP

Sizing chart

All measurements are shown in cm


Below knee

	S	M	L	XL
Ankle (cB)	19-23	23-26	26-29	29-32
Calf (cC)	31-36	34-39	37-42	40-45
EMIS TIP	SMA	MED	LAR	LAR

Length


	Short	Long
Below knee (lD)	33-37	38-44
EMIS TIP	NOR	LON

Below knee



Styles and colours

● Closed toe



Beige

	EMIS TIP	
ULCER X Liner	CLO	●




Scan here for this product's ordering codes

ULCER X




Product codes

Below knee



		S		M		L		XL	
		STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS
ULCER X Liners	Short	81080	-	81082	-	81084	-	81086	-
	Long	81081	-	81083	-	81085	-	81087	-
ULCER X kit	Short	81064	81072	81066	81074	81068	81076	81070	81078
	Long	81065	81073	81067	81075	81069	81077	81071	81079



Sensinnov
grip top band
on thigh



Thigh with grip top
in *Marine*



Below knee in
Dune



Tights (panty) in
Black

Material composition

Below knee

72% Polyamide, 28% Elastane

Thigh

64% Polyamide,
28% Elastane, 8% Polyester

Tights

72% Polyamide,
27% Elastane, 1% Cotton
All garments latex free



Scan here
for this product's
measuring form

Essential SEMITRSPARENT

Fashionable, comfortable, semi-transparent compression garments for everyday wear, with the advantage of easy donning and doffing.

Compression class

10-15 mmHg	15-20 mmHg	20-36 mmHg
AFNOR 1	AFNOR 2	AFNOR 3

Product benefits

Easy donning & doffing	Elegant	Soft

Suitable for

- Venous Insufficiency
- Mild lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Lipoedema (lipalgia)
- Mixed disease

EMIS search tips

Product name

Essential SEMITRSPARENT

SIG ESS SEMI
EMIS TIP

Compression class (mmHg level)

Class 1 10 EMIS TIP Class 2 15 EMIS TIP Class 3 20 EMIS TIP

Sizing chart

All measurements are shown in cm

Below knee

	XS*	S		M		L		XL		XXL*
Ankle (cB)	17-19	19-22		22-24		24-26		26-29		29-32
	STD	SLIM*	STD	PLUS*	SLIM*	STD	PLUS*	SLIM*	STD	PLUS*
Calf (cC)	27-35	22-30	30-38	38-46	25-33	33-41	41-49	29-37	37-45	45-53
EMIS TIP	SMA	SMA SLI	SMA	SMA PLU	MED SLI	MED	MED PLU	LAR SLI	LAR	LAR PLU


Thigh with grip top / tights (panty)

	XS*	S		M		L		XL		XXL*
Ankle (cB)	17-19	19-22		22-24		24-26		26-29		29-32
	STD	SLIM*	STD	PLUS*	SLIM*	STD	PLUS*	SLIM*	STD	PLUS*
Calf (cC)	27-35	26-34	30-38	34-42	29-37	33-41	37-45	33-41	37-45	41-49
Thigh (cG)	42-55	35-45	45-58	58-68	39-49	49-62	62-72	44-54	54-67	67-77
Hip (cH) (For tights only)	-	-	≤100	≤108‡	-	≤107	≤118‡	-	≤114	≤128‡
EMIS TIP	SMA	SMA SLI	SMA	SMA PLU	MED SLI	MED	MED PLU	LAR SLI	LAR	LAR PLU

‡ Tights plus size only available in Compression Class 2 * These products and sizes are only available in Compression Class 2 & 3

Length

	Normal	Long
Below knee (ℓD)	≤ 40	≥ 40
Thigh (ℓG)	≤ 75	≥ 75
Tights (panty) (ℓG)	≤ 75	≥ 75
Shoe size	36/40 - 3.5/6.5	37/41 - 4/7.5
EMIS TIP	NOR	LON



Styles and colours

	EMIS TIP	Class 1		Class 2				Class 3			
		Black	Dune	Black	Dune	Grey ≈	Marine ≈	Natural	Spice ≈	Black	Dune
Below knee	BEL	●	●	●	●	●	●	●	●	●	●
Thigh with grip top	THI	●	●	●	●	●	●	●	●	●	●
Tights (panty)	PAN	●	●	●	●	●	●	●	●	●	●
EMIS TIP		BLA	DUN	BLA	DUN	GRE	MAR	NAT	SPI	BLA	DUN

≈ Grey, Marine and Spice colours only available in standard size



Scan here
for this
product's
ordering
codes

Essential SEMITRSPARENT



Product codes

Below knee




			Class 1		Class 2					Class 3		
			Black	Dune	Black	Dune	Grey	Marine	Natural	Spice	Black	Dune
			●	●	●	●	●	●	●	●	●	●
XS	STD	Normal	-	-	312610	312624	313319	312695	312597	312637	313368	313380
		Long	-	-	312617	312631	313289	312701	312604	312643	313374	313386
S	SLIM	Normal	-	-	312776	312785	-	-	312767	-	313404	313412
		Long	-	-	312780	312789	-	-	312771	-	313408	313416
	STD	Normal	313320	313328	312611	312625	313284	312696	312598	312638	313369	313381
		Long	313324	313332	312618	312632	313290	312702	312605	312644	313375	313387
	PLUS	Normal	-	-	312803	312813	-	-	312793	-	313421	313429
		Long	-	-	312808	312818	-	-	312798	-	313425	313433
M	SLIM	Normal	-	-	312777	312786	-	-	312768	-	313405	313413
		Long	-	-	312781	312790	-	-	312772	-	313409	313417
	STD	Normal	313321	313329	312612	312626	313285	312697	312599	312639	313370	313382
		Long	313325	313333	312619	312633	313291	312703	312606	312645	313376	313388
	PLUS	Normal	-	-	312804	312814	-	-	312794	-	313422	313430
		Long	-	-	312809	312820	-	-	312799	-	313426	313434
L	SLIM	Normal	-	-	312778	312787	-	-	312769	-	313406	313414
		Long	-	-	312783	312791	-	-	312773	-	313410	313419
	STD	Normal	313322	313330	312613	312627	313286	312698	312600	312640	313371	313383
		Long	313326	313334	312620	312634	313292	312704	312607	312646	313377	313389
	PLUS	Normal	-	-	312805	312815	-	-	312795	-	313423	313431
		Long	-	-	312810	312821	-	-	312800	-	313427	313435
XL	SLIM	Normal	-	-	312779	312788	-	-	312770	-	313407	313415
		Long	-	-	312784	312792	-	-	312774	-	313411	313420
	STD	Normal	313323	313331	312614	312629	313287	312699	312602	312641	313372	313384
		Long	313327	313335	312684	312635	313293	312705	312608	312647	313378	313390
	PLUS	Normal	-	-	312806	312816	-	-	312796	-	313424	313432
		Long	-	-	312811	312822	-	-	312801	-	313428	313436
XXL	STD	Normal	-	-	312615	312630	313288	312700	312603	312642	313373	313385
		Long	-	-	312685	312636	313294	312706	312616	312648	313379	313391
	PLUS	Normal	-	-	312807	312817	-	-	312797	-	-	-
		Long	-	-	312812	312823	-	-	312802	-	-	-

●●●●● Closed toe

Product codes

Thigh with grip top




			Class 1		Class 2					Class 3		
			Black	Dune	Black	Dune	Grey	Marine	Natural	Spice	Black	Dune
			●	●	●	●	●	●	●	●	●	●
XS	STD	Normal	-	-	312846	312871	313295	312859	312834	312904	313437	313449
		Long	-	-	312853	312877	313301	312865	312840	312910	313443	313455
S	SLIM	Normal	-	-	313018	313026	-	-	313034	-	313481	313473
		Long	-	-	313022	313030	-	-	313038	-	313485	313477
	STD	Normal	313336	313344	312848	312872	313296	312860	312835	312905	313438	313450
		Long	313340	313348	312854	312878	313302	312866	312841	312911	313444	313456
	PLUS	Normal	-	-	313052	313062	-	-	313042	-	313497	313489
		Long	-	-	313057	313067	-	-	313047	-	313501	313493
M	SLIM	Normal	-	-	313019	313027	-	-	313035	-	313482	313474
		Long	-	-	313023	313031	-	-	313039	-	313486	313478
	STD	Normal	313337	313345	312849	312873	313297	312861	312836	312906	313439	313451
		Long	313341	313349	312855	312879	313303	312867	312842	312912	313445	313457
	PLUS	Normal	-	-	313053	313063	-	-	313043	-	313498	313490
		Long	-	-	313058	313068	-	-	313048	-	313502	313494
L	SLIM	Normal	-	-	313020	313028	-	-	313036	-	313483	313475
		Long	-	-	313024	313032	-	-	313040	-	313487	313479
	STD	Normal	313338	313346	312850	312874	313298	312862	312837	312907	313440	313452
		Long	313342	313350	312856	312880	313304	312868	312843	312913	313446	313458
	PLUS	Normal	-	-	313054	313064	-	-	313044	-	313499	313491
		Long	-	-	313059	313069	-	-	313049	-	313503	313495
XL	SLIM	Normal	-	-	313021	313029	-	-	313037	-	313484	313476
		Long	-	-	313025	313033	-	-	313041	-	313488	313480
	STD	Normal	313339	313347	312851	312875	313299	312863	312838	312908	313441	313453
		Long	313343	313351	312857	312881	313305	312869	312844	312914	313447	313459
	PLUS	Normal	-	-	313055	313065	-	-	313045	-	313500	313492
		Long	-	-	313060	313070	-	-	313050	-	313504	313496
XXL	STD	Normal	-	-	312852	312876	313300	312864	312839	312909	313442	313454
		Long	-	-	312858	312882	313306	312870	312845	312915	313448	313460
	PLUS	Normal	-	-	313056	313066	-	-	313046	-	-	-
		Long	-	-	313061	313071	-	-	313051	-	-	-

●●●●● Closed toe

Product codes

Tights (panty)




			Class 1		Class 2				Class 3			
			Black	Dune	Black	Dune	Grey	Marine	Natural	Spice	Black	Dune
			●	●	●	●	●	●	●	●	●	●
S	STD	Normal	313352	313360	313104	313135	313307	313125	313094	313145	313506	313514
		Long	313356	313364	313109	313140	313312	313130	313099	313150	313510	313518
	PLUS	Normal	-	-	313155	313176	-	-	313165	-	-	-
		Long	-	-	313160	313181	-	-	313171	-	-	-
M	STD	Normal	313353	313361	313105	313136	313308	313126	313095	313146	313507	313515
		Long	313357	313365	313110	313141	313313	313131	313100	313151	313511	313519
	PLUS	Normal	-	-	313156	313177	-	-	313166	-	-	-
		Long	-	-	313161	313182	-	-	313172	-	-	-
L	STD	Normal	313354	313362	313106	313137	313309	313127	313096	313147	313508	313516
		Long	313358	313366	313111	313142	313314	313132	313101	313152	313512	313520
	PLUS	Normal	-	-	313157	313178	-	-	313168	-	-	-
		Long	-	-	313162	313183	-	-	313173	-	-	-
XL	STD	Normal	313355	313363	313107	313138	313310	313128	313097	313148	313509	313517
		Long	313359	3131367	313112	313143	313315	313133	313102	313153	313513	313521
	PLUS	Normal	-	-	313158	313179	-	-	313169	-	-	-
		Long	-	-	313163	313184	-	-	313174	-	-	-
XXL	STD	Normal	-	-	313108	313139	313311	313129	313098	313149	-	-
		Long	-	-	313113	313144	313316	313134	313103	313154	-	-
	PLUS	Normal	-	-	313159	313180	-	-	313170	-	-	-
		Long	-	-	313164	313185	-	-	313175	-	-	-

●●●●● Closed toe



Open toe
on below knee
and thigh



SensinnoV
grip top band
on thigh



Below knee
in Beige
Rose



Thigh with grip top
in Black

Material composition

Below knee

78% Polyamide, 22% Elastane

Thigh

72% Polyamide,
24% Elastane, 4% Polyester
All garments latex free



Scan here
for this product's
measuring form

Style OPAQUE

A soft, fully opaque fabric for discreet, effective management of mild to moderate conditions. Ideal for long-term wear in cooler climates or professional settings.

Compression class

15-20
mmHg

AFNOR
2

Product benefits



Opaque

Easy
donning &
doffing

Soft

Suitable for

- Venous Insufficiency
- Mild lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Lipoedema (lipalgia)
- Mixed disease

EMIS search tips

Product name

Style OPAQUE

SIG STY OPA

EMIS TIP

Sizing chart

All measurements are shown in cm

Below knee

	XS	S		M		L		XL		XXL					
Ankle (cB)	17-19	19-22		22-24		24-26		26-29		29-32					
	STD	SLIM	STD	PLUS	SLIM	STD	PLUS	SLIM	STD	PLUS	SLIM	STD	PLUS		
Calf (cC)	27-35	22-30	30-38	38-46	25-33	33-41	41-49	29-37	37-45	45-53	32-40	40-48	48-56	43-51	51-59
EMIS TIP	SMA	SMA SLI	SMA	SMA PLU	MED SLI	MED	MED PLU	LAR SLI	LAR	LAR PLU	LAR SLI	LAR	LAR PLU	LAR	LAR PLU


Thigh with grip top

	XS	S		M		L		XL		XXL					
Ankle (cB)	17-19	19-22		22-24		24-26		26-29		29-32					
	STD	SLIM	STD	PLUS	SLIM	STD	PLUS	SLIM	STD	PLUS	SLIM	STD	PLUS		
Calf (cC)	27-35	26-34	30-38	34-42	29-37	33-41	37-45	33-41	37-45	41-49	36-44	40-48	44-52	43-51	47-55
Thigh (cG)	42-55	35-45	45-58	58-68	39-49	49-62	62-72	44-54	54-67	67-77	46-56	56-69	69-79	62-75	75-85
EMIS TIP	SMA	SMA SLI	SMA	SMA PLU	MED SLI	MED	MED PLU	LAR SLI	LAR	LAR PLU	LAR SLI	LAR	LAR PLU	LAR	LAR PLU

Length

	Normal	Long
Below knee (lD)	≤ 40	≥ 40
Thigh (lG)	≤ 75	≥ 75
Shoe size	36/40 - 3.5/6.5	37/41 - 4/7.5
EMIS TIP	NOR	LON

Below knee **Thigh**



Styles and colours

		Class 2		
		Beige Rose*	Black	Marine**
Below knee	BEL	● ●	● ●	●
Thigh with grip top	THI	● ●	● ●	●
EMIS TIP		BEI	BLA	MAR

* Beige Rose is not available in plus sizes.

** Marine and open toe products are only available in standard sizes.



Scan here
for this
product's
ordering
codes

Style OPAQUE



Product codes

Below knee




			Class 2				
			Beige Rose		Black		Marine
			●	○	●	○	●
XS	STD	Normal	301460	301460-OT	301448	301448-OT	301472
		Long	301466	301466-OT	301454	301454-OT	301478
S	SLIM	Normal	301492	-	301484	-	-
		Long	301496	-	301488	-	-
	STD	Normal	301461	301461-OT	301449	301449-OT	301473
		Long	301467	301467-OT	301455	301455-OT	301479
	PLUS	Normal	-	-	301500	-	-
		Long	-	-	301505	-	-
M	SLIM	Normal	301493	-	301485	-	-
		Long	301497	-	301489	-	-
	STD	Normal	301462	301462-OT	301450	301450-OT	301474
		Long	301468	301468-OT	301456	301456-OT	301480
	PLUS	Normal	-	-	301501	-	-
		Long	-	-	301506	-	-
L	SLIM	Normal	301494	-	301486	-	-
		Long	301498	-	301490	-	-
	STD	Normal	301463	301463-OT	301451	301451-OT	301475
		Long	301469	301469-OT	301457	301457-OT	301481
	PLUS	Normal	-	-	301502	-	-
		Long	-	-	301507	-	-
XL	SLIM	Normal	301495	-	301487	-	-
		Long	301499	-	301491	-	-
	STD	Normal	301464	301464-OT	301452	301452-OT	301476
		Long	301470	301470-OT	301458	301458-OT	301482
	PLUS	Normal	-	-	301503	-	-
		Long	-	-	301508	-	-
XXL	STD	Normal	301465	301465-OT	301453	301453-OT	301477
		Long	301471	301471-OT	301459	301459-OT	301483
	PLUS	Normal	-	-	301504	-	-
		Long	-	-	301509	-	-

●●● Closed toe ○○ Open toe

Product codes

Thigh with grip top



			Class 2				
			Beige Rose		Black		Marine
			●	○	●	○	●
XS	STD	Normal	301522	301522-OT	301510	301510-OT	301534
		Long	301528	301528-OT	301516	301516-OT	301540
S	SLIM	Normal	301554	-	301546	-	-
		Long	301558	-	301550	-	-
	STD	Normal	301523	301523-OT	301511	301511-OT	301535
		Long	301529	301529-OT	301517	301517-OT	301541
	PLUS	Normal	-	-	301562	-	-
		Long	-	-	301567	-	-
M	SLIM	Normal	301555	-	301547	-	-
		Long	301559	-	301551	-	-
	STD	Normal	301524	301524-OT	301512	301512-OT	301536
		Long	301530	301530-OT	301518	301518-OT	301542
	PLUS	Normal	-	-	301563	-	-
		Long	-	-	301568	-	-
L	SLIM	Normal	301556	-	301548	-	-
		Long	301560	-	301552	-	-
	STD	Normal	301525	301525-OT	301513	301513-OT	301537
		Long	301531	301531-OT	301519	301519-OT	301543
	PLUS	Normal	-	-	301564	-	-
		Long	-	-	301569	-	-
XL	SLIM	Normal	301557	-	301549	-	-
		Long	301561	-	301553	-	-
	STD	Normal	301526	301526-OT	301514	301514-OT	301538
		Long	301532	301532-OT	301520	301520-OT	301544
	PLUS	Normal	-	-	301565	-	-
		Long	-	-	301570	-	-
XXL	STD	Normal	301527	301527-OT	301515	301515-OT	301539
		Long	301533	301533-OT	301521	301521-OT	301545
	PLUS	Normal	-	-	301566	-	-
		Long	-	-	301571	-	-

●●● Closed toe ○○ Open toe

Available in MTM

Offers tailored sizing, colour options, and style choices, ensuring the right product with the right fit.

Find out more in the Made to Measure section

Material composition

Below knee

61% Polyamide, 39% Elastane

Thigh

63% Polyamide, 37% Elastane

Tights

59% Polyamide, 41% Elastane

All garments latex free



Scan here for this product's measuring form

Style SEMITRSPARENT

A breathable, semi-transparent design, providing an effective therapy for venous conditions with a discreet, elegant finish to support patient comfort.

Compression class

18-21 mmHg

RAL 1

23-32 mmHg

RAL 2


Product benefits



Semi-transparent, Durable, Breathable, Made to Measure

Suitable for

- Venous insufficiency
- Mild chronic oedema
- Mixed disease
- Lipoedema (lipalgia)
- Prevention of leg ulcers and oedema
- Pre & post vein treatment with minimal oedema



Sensinnov laced grip top on thigh



Knobbed grip top band on thigh



Thigh with lace grip top in Skin



Below knee in Black with Open Toe



Tights (panty) in Skin

EMIS search tips

Product name

Style SEMITRSPARENT

SIG STY SEMI
EMIS TIP

Compression class (mmHg level)

Class 1 18
EMIS TIP

Class 2 23
EMIS TIP

Sizing chart


All measurements are shown in cm

Below knee, thigh and tights (panty)

	XXS		XS		S		M		L		XL	
Ankle (cB)	18-20		20-22		22-24		24-26		26-29		29-32	
	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS
Calf (cC)	25-32	29-36	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49	42-49	43-51
Thigh (cG) with grip top	42-52	47-56	47-56	52-61	52-61	56-65	56-65	61-71	61-71	66-76	66-76	73-85
Thigh (cG) (For tights only)	42-54	47-59	47-59	52-64	52-64	56-68	56-68	61-74	61-74	66-79	66-79	76-85
Hip (cH)	75-90	80-95	80-95	85-100	85-100	95-110	95-110	105-120	105-120	115-130	115-130	125-140
Waist (cT)	60-75	65-80	65-80	75-90	75-90	85-100	85-100	95-110	95-110	105-120	105-120	115-130
EMIS TIP	SMA	SMA PLU	SMA	SMA PLU	SMA	SMA PLU	MED	MED PLU	LAR	LAR PLU	LAR	LAR PLU

Length

	Short	Normal	Long
Below knee (ℓD)	-	35-40	41-45
Thigh (ℓG)	56-64	65-74	75-84
Tights (panty) (ℓG)	-	65-74	75-84
EMIS TIP	SHO	NOR	LON



Styles and colours

- Closed toe CLO
 - Open toe OPE
- EMIS TIP

	EMIS TIP	Class 1		Class 2	
		Black	Skin	Black	Skin
Below knee	BEL	● ○	○ ○	● ○	○ ○
Thigh with lace grip top	THI GRI	● ○	○ ○	● ○	○ ○
Thigh with knobbed grip top	THI GRI	● ○	○ ○	● ○	○ ○
Tights (panty)	PAN	● ○	○ ○	● ○	○ ○
EMIS TIP		BLK	SKI	BLK	SKI



Scan here for this product's ordering codes

Style SEMITRSPARENT



Product codes

Below knee




			Class 1				Class 2			
			Black		Skin		Black		Skin	
			●	○	●	○	●	○	●	○
XXS	STD	Normal	54172	54170	54168	54166	54264	54262	54260	54258
		Long	54173	54171	54169	54167	54265	54263	54261	54259
	PLUS	Normal	66063	66059	66055	66051	55166	55160	55154	55148
		Long	66064	66060	66056	66052	55167	55161	55155	55149
XS	STD	Normal	35042	34898	35018	34874	36482	36338	36458	36314
		Long	35043	34899	35019	34875	36483	36339	36459	36315
	PLUS	Normal	35054	34910	35030	34886	36494	36350	36470	36326
		Long	35055	34911	35031	34887	36495	36351	36471	36327
S	STD	Normal	35045	34901	35021	34877	36485	36341	36461	36317
		Long	35046	34902	35022	34878	36486	36342	36462	36318
	PLUS	Normal	35057	34913	35033	34889	36497	36353	36473	36329
		Long	35058	34914	35034	34890	36498	36354	36474	36330
M	STD	Normal	35048	34904	35024	34880	36488	36344	36464	36320
		Long	35049	34905	35025	34881	36489	36345	36465	36321
	PLUS	Normal	35060	34916	35036	34892	36500	36356	36476	36332
		Long	35061	34917	35037	34893	36501	36357	36477	36333
L	STD	Normal	35051	34907	35027	34883	36491	36347	36467	36323
		Long	35052	34908	35028	34884	36492	36348	36468	36324
	PLUS	Normal	35063	34919	35039	34895	36503	36359	36479	36335
		Long	35064	34920	35040	34896	36504	36360	36480	36336
XL	STD	Normal	55842	55840	55838	55836	55164	55158	55152	55146
		Long	55843	55841	55839	55837	55165	55159	55153	55147
	PLUS	Normal	66065	66061	66057	66053	55168	55162	55156	55150
		Long	66066	66062	66058	66054	55169	55163	55157	55151

● Closed toe ○ Open toe

Product codes

Thigh with lace grip top




			Class 1				Class 2			
			Black		Skin		Black		Skin	
			●	○	●	○	●	○	●	○
XXS	STD	Short	58652	58616	58580	58544	57430	57394	57358	57322
		Normal	58653	58617	58581	58545	57431	57395	57359	57323
		Long	58654	58618	58582	58546	57432	57396	57360	57324
	PLUS	Short	58670	58634	58598	58562	57448	57412	57376	57340
		Normal	58671	58635	58599	58563	57449	57413	57377	57341
		Long	58672	58636	58600	58564	57450	57414	57378	57342
XS	STD	Short	58655	58619	58583	58547	57433	57397	57361	57325
		Normal	58656	58620	58584	58548	57434	57398	57362	57326
		Long	58657	58621	58585	58549	57435	57399	57363	57327
	PLUS	Short	58673	58637	58601	58565	57451	57415	57379	57343
		Normal	58674	58638	58602	58566	57452	57416	57380	57344
		Long	58675	58639	58603	58567	57453	57417	57381	57345
S	STD	Short	58658	58622	58586	58550	57436	57400	57364	57328
		Normal	58659	58623	58587	58551	57437	57401	57365	57329
		Long	58660	58624	58588	58552	57438	57402	57366	57330
	PLUS	Short	58676	58640	58604	58568	57454	57418	57382	57346
		Normal	58677	58641	58605	58569	57455	57419	57383	57347
		Long	58678	58642	58606	58570	57456	57420	57384	57348
M	STD	Short	58661	58625	58589	58553	57439	57403	57367	57331
		Normal	58662	58626	58590	58554	57440	57404	57368	57332
		Long	58663	58627	58591	58555	57441	57405	57369	57333
	PLUS	Short	58679	58643	58607	58571	57457	57421	57385	57349
		Normal	58680	58644	58608	58572	57458	57422	57386	57350
		Long	58681	58645	58609	58573	57459	57423	57387	57351
L	STD	Short	58664	58628	58592	58556	57442	57406	57370	57334
		Normal	58665	58629	58593	58557	57443	57407	57371	57335
		Long	58666	58630	58594	58558	57444	57408	57372	57336
	PLUS	Short	58682	58646	58610	58574	57460	57424	57388	57352
		Normal	58683	58647	58611	58575	57461	57425	57389	57353
		Long	58684	58648	58612	58576	57462	57426	57390	57354
XL	STD	Short	58667	58631	58595	58559	57445	57409	57373	57337
		Normal	58668	58632	58596	58560	57446	57410	57374	57338
		Long	58669	58633	58597	58561	57447	57411	57375	57339
	PLUS	Short	58685	58649	58613	58577	57463	57427	57391	57355
		Normal	58686	58650	58614	58578	57464	57428	57392	57356
		Long	58687	58651	58615	58579	57465	57429	57393	57357

● Closed toe ○ Open toe

Product codes

Thigh with knobbed grip top




			Class 1				Class 2			
			Black		Skin		Black		Skin	
			●	○	●	○	●	○	●	○
XXS	STD	Short	54226	54223	54220	54217	54275	54272	54269	54266
		Normal	54227	54224	54221	54218	54276	54273	54270	54267
		Long	54228	54225	54222	54219	54277	54274	54271	54268
	PLUS	Short	66305	66296	66287	66278	55493	55484	55475	55466
		Normal	66306	66297	66288	66279	55494	55485	55476	55467
		Long	66307	66298	66289	66280	55495	55486	55477	55468
XS	STD	Short	35617	35473	35593	35449	37057	36913	37033	36889
		Normal	35618	35474	35594	35450	37058	36914	37034	36890
		Long	35619	35475	35595	35451	37059	36915	37035	36891
	PLUS	Short	35629	35485	35605	35461	37069	36925	37045	36901
		Normal	35630	35486	35606	35462	37070	36926	37046	36902
		Long	35631	35487	35607	35463	37071	36927	37047	36903
S	STD	Short	35620	35476	35596	35452	37060	36916	37036	36892
		Normal	35621	35477	35597	35453	37061	36917	37037	36893
		Long	35622	35478	35598	35454	37062	36918	37038	36894
	PLUS	Short	35632	35488	35608	35464	37072	36928	37048	36904
		Normal	35633	35489	35609	35465	37073	36929	37049	36905
		Long	35634	35490	35610	35466	37074	36930	37050	36906
M	STD	Short	35623	35479	35599	35455	37063	36919	37039	36895
		Normal	35624	35480	35600	35456	37064	36920	37040	36896
		Long	35625	35481	35601	35457	37065	36921	37041	36897
	PLUS	Short	35635	35491	35611	35467	37075	36931	37051	36907
		Normal	35636	35492	35612	35468	37076	36932	37052	36908
		Long	35637	35493	35613	35469	37077	36933	37053	36909
L	STD	Short	35626	35482	35602	35458	37066	36922	37042	36898
		Normal	35627	35483	35603	35459	37067	36923	37043	36899
		Long	35628	35484	35604	35460	37068	36924	37044	36900
	PLUS	Short	35638	35494	35614	35470	37078	36934	37054	36910
		Normal	35639	35495	35615	35471	37079	36935	37055	36911
		Long	35640	35496	35616	35472	37080	36936	37056	36912
XL	STD	Short	66302	66293	66284	66275	55490	55481	55472	55463
		Normal	66303	66294	66285	66276	55491	55482	55473	55464
		Long	66304	66295	66286	66277	55492	55483	55474	55465
	PLUS	Short	66308	66299	66290	66281	55496	55487	55478	55469
		Normal	66309	66300	66291	66282	55497	55488	55479	55470
		Long	66310	66301	66292	66283	55498	55489	55480	55471

●● Closed toe ○○ Open toe

Product codes

Tights (panty)



			Class 1				Class 2			
			Black		Skin		Black		Skin	
			●	○	●	○	●	○	●	○
XXS	STD	Normal	54248	54246	54244	54242	54296	54294	54292	54290
		Long	54249	54247	54245	54243	54297	54295	54293	54291
	PLUS	Normal	66747	66743	66739	66735	55648	55642	55636	55630
		Long	66748	66744	66740	66736	55649	55643	55637	55631
XS	STD	Normal	35906	35762	35882	35738	37346	37202	37322	37178
		Long	35907	35763	35883	35739	37347	37203	37323	37179
	PLUS	Normal	35918	35774	35894	35750	37358	37214	37334	37190
		Long	35919	35775	35895	35751	37359	37215	37335	37191
S	STD	Normal	35909	35765	35885	35741	37349	37205	37325	37181
		Long	35910	35766	35886	35742	37350	37206	37326	37182
	PLUS	Normal	35921	35777	35897	35753	37361	37217	37337	37193
		Long	35922	35778	35898	35754	37362	37218	37338	37194
M	STD	Normal	35912	35768	35888	35744	37352	37208	37328	37184
		Long	35913	35769	35889	35745	37353	37209	37329	37185
	PLUS	Normal	35924	35780	35900	35756	37364	37220	37340	37196
		Long	35925	35781	35901	35757	37365	37221	37341	37197
L	STD	Normal	35915	35771	35891	35747	37355	37211	37331	37187
		Long	35916	35772	35892	35748	37356	37212	37332	37188
	PLUS	Normal	35927	35783	35903	35759	37367	37223	37343	37199
		Long	35928	35784	35904	35760	37368	37224	37344	37200
XL	STD	Normal	55862	55860	55858	55856	55646	55640	55634	55628
		Long	55863	55861	55859	55857	55647	55641	55635	55629
	PLUS	Normal	66749	66745	66741	66737	55650	55644	55638	55632
		Long	66750	66746	66742	66738	55651	55645	55639	55633

●● Closed toe ○○ Open toe



Open toe available on some styles



Sensinnov grip top band on thigh



Thigh with grip top in Beige Four (140)



Below knee in Beige One (110)



Tights (panty) in Black

Material composition

Below knee

61% Polyamide, 38% Elastane, 1% Polyester

Thigh

58% Polyamide, 42% Elastane

Tights

62% Polyamide, 37% Elastane, 1% Polyester

All garments latex free



Scan here for this product's measuring form

Style TRANSPARENT

Blends a fashionable sheer look with long lasting softness for professional and everyday wear.

Compression class

15-20 mmHg

AFNOR 2

Product benefits



Nude



Microfibre



Reinforced sole

Suitable for

- Venous insufficiency
- Mild lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Lipoedema (lipalgia)
- Mixed disease

EMIS search tips

Product name
Style TRANSPARENT

SIG STY TRANS
EMIS TIP

Compression class (mmHg level)

Class 2 15
EMIS TIP

Sizing chart

All measurements are shown in cm

Below knee

	XS	S		M		L		XL		XXL
Ankle (cB)	17-19	19-22		22-24		24-26		26-29		29-32
	STD	SLIM	STD	PLUS	SLIM	STD	PLUS	SLIM	STD	PLUS
Calf (cC)	27-35	22-30	30-38	38-46	25-33	33-41	41-49	29-37	37-45	45-53
	32-40	40-48	48-56	43-51	51-59					
EMIS TIP	SMA	SMA SLI	SMA	SMA PLU	MED SLI	MED	MED PLU	LAR SLI	LAR	LAR PLU


Thigh with grip top / tights (panty)

	XS	S		M		L		XL		XXL
Ankle (cB)	17-19	19-22		22-24		24-26		26-29		29-32
	STD	SLIM	STD	PLUS	SLIM	STD	PLUS	SLIM	STD	PLUS
Calf (cC)	27-35	26-34	30-38	34-42	29-37	33-41	37-45	33-41	37-45	41-49
	36-44	40-48	44-52	43-51	47-55					
Thigh (cG)	42-55	35-45	45-58	58-68	39-49	49-62	62-72	44-54	54-67	67-77
	46-56	56-69	69-79	62-75	75-85					
Hip (cH) (For tights only)	-	-	≤108	108≥	-	≤118	118≥	-	≤128	128≥
	-	-	≤138	138≥	≤138	138≥				
EMIS TIP	SMA	SMA SLI	SMA	SMA PLU	MED SLI	MED	MED PLU	LAR SLI	LAR	LAR PLU


Length

	Normal	Long
Below knee (ℓD)	≤ 40	≥ 40
Thigh (ℓG)	≤ 75	≥ 75
Tights (panty) (ℓG)	≤ 75	≥ 75
Shoe size	36/40 - 3.5/6.5	37/41 - 4/7.5
EMIS TIP	NOR	LON


Below knee



Thigh



Tights (panty)



Styles and colours

		Class 2						
		Beige One (110)*	Beige Two (120)*	Beige Three (130)*	Beige Four (140)	Beige Five (150)	Black	Midnight Blue
●●● Closed toe	CLO							
○ Open toe	OPE							
Below knee	BEL	● ○	● ○	● ○	● ○	●	●	●
Thigh with grip top	THI	● ○	● ○	● ○	● ○	●	●	●
Tights (panty)	PAN	● ○	● ○	● ○	● ○	●	●	●
EMIS TIP		ONE	TWO	THREE	FOUR	FIVE	BLA	MID

* Slim and plus size only available in Beige One, Two and Three



Scan here for this product's ordering codes

Style TRANSPARENT



Product codes

Below knee




			Class 2										
			Beige One (110)		Beige Two (120)		Beige Three (130)		Beige Four (140)		Beige Five (150)	Black	Midnight Blue
			●	○	●	○	●	○	●	○	●	●	●
XS	STD	Normal	301612	301738	301624	301750	301636	307607	301648	301762	301660	301672	313614
		Long	301618	301744	301630	301756	301642	307613	301654	301768	301666	301678	313620
S	SLIM	Normal	301684	-	301692	-	307598	-	-	-	-	301700	-
		Long	301688	-	301696	-	307602	-	-	-	-	301704	-
	STD	Normal	301613	301739	301625	301751	301637	307608	301649	301763	301661	301673	313615
		Long	301619	301745	301631	301757	301643	307614	301655	301769	301667	301679	313621
	PLUS	Normal	301708	-	301718	-	307559	-	-	-	-	301728	-
		Long	301713	-	301723	-	307564	-	-	-	-	301733	-
M	SLIM	Normal	301685	-	301693	-	307599	-	-	-	-	301701	-
		Long	301689	-	301697	-	307604	-	-	-	-	301705	-
	STD	Normal	301614	301740	301626	301752	301638	307609	301650	301764	301662	301674	313616
		Long	301620	301746	301632	301758	301644	307615	301656	301770	301668	301680	313622
	PLUS	Normal	301709	-	301719	-	307560	-	-	-	-	301729	-
		Long	301714	-	301724	-	307565	-	-	-	-	301734	-
L	SLIM	Normal	301686	-	301694	-	307600	-	-	-	-	301702	-
		Long	301690	-	301698	-	307605	-	-	-	-	301706	-
	STD	Normal	301615	301741	301627	301753	301639	307610	301651	301765	301663	301675	313617
		Long	301621	301747	301633	301759	301645	307616	301657	301771	301669	301681	313623
	PLUS	Normal	301710	-	301720	-	307561	-	-	-	-	301730	-
		Long	301715	-	301725	-	307566	-	-	-	-	301735	-
XL	SLIM	Normal	301687	-	301695	-	307601	-	-	-	-	301703	-
		Long	301691	-	301699	-	307606	-	-	-	-	301707	-
	STD	Normal	301616	301742	301628	301754	301640	307611	301652	301766	301664	301676	313618
		Long	301622	301748	301634	301760	301646	307617	301658	301772	301670	301682	313624
	PLUS	Normal	301711	-	301721	-	307562	-	-	-	-	301731	-
		Long	301716	-	301726	-	307567	-	-	-	-	301736	-
XXL	STD	Normal	301617	301743	301629	301755	301641	307612	301653	301767	301665	301677	313619
		Long	301623	301749	301635	301761	301647	307618	301659	301773	301671	301683	313625
	PLUS	Normal	301712	-	301722	-	307563	-	-	-	-	301732	-
		Long	301717	-	301727	-	307568	-	-	-	-	301737	-

●●● Closed toe ○ Open toe

Product codes

Thigh with grip top




			Class 2										
			Beige One (110)		Beige Two (120)		Beige Three (130)		Beige Four (140)		Beige Five (150)	Black	Midnight Blue
			●	○	●	○	●	○	●	○	●	●	●
XS	STD	Normal	301774	301900	301786	301912	301798	307619	301810	301924	301822	301834	313627
		Long	301780	301906	301792	301918	301804	307625	301816	301930	301828	301840	313633
S	SLIM	Normal	301846	-	301854	-	307569	-	-	-	-	301862	-
		Long	301850	-	301858	-	307573	-	-	-	-	301866	-
	STD	Normal	301775	301901	301787	301913	301799	307620	301811	301925	301823	301835	313628
		Long	301781	301907	301793	301919	301805	307626	301817	301931	301829	301841	313634
	PLUS	Normal	301870	-	301880	-	307577	-	-	-	-	301890	-
		Long	301875	-	301885	-	307582	-	-	-	-	301895	-
M	SLIM	Normal	301847	-	301855	-	307570	-	-	-	-	301863	-
		Long	301851	-	301859	-	307574	-	-	-	-	301867	-
	STD	Normal	301776	301902	301788	301914	301800	307621	301812	301926	301824	301836	313629
		Long	301782	301908	301794	301920	301806	307627	301818	301932	301830	301842	313635
	PLUS	Normal	301871	-	301881	-	307578	-	-	-	-	301891	-
		Long	301876	-	301886	-	307583	-	-	-	-	301896	-
L	SLIM	Normal	301848	-	301856	-	307571	-	-	-	-	301864	-
		Long	301852	-	301860	-	307575	-	-	-	-	301868	-
	STD	Normal	301777	301903	301789	301915	301801	307622	301813	301927	301825	301837	313630
		Long	301783	301909	301795	301921	301807	307628	301819	301933	301831	301843	313636
	PLUS	Normal	301872	-	301882	-	307579	-	-	-	-	301892	-
		Long	301877	-	301887	-	307584	-	-	-	-	301897	-
XL	SLIM	Normal	301849	-	301857	-	307572	-	-	-	-	301865	-
		Long	301853	-	301861	-	307576	-	-	-	-	301869	-
	STD	Normal	301778	301904	301790	301916	301802	307623	301814	301928	301826	301838	313631
		Long	301784	301910	301796	301922	301808	307630	301820	301934	301832	301844	313637
	PLUS	Normal	301873	-	301883	-	307580	-	-	-	-	301893	-
		Long	301878	-	301888	-	307585	-	-	-	-	301898	-
XXL	STD	Normal	301779	301905	301791	301917	301803	307624	301815	301929	301827	301839	313632
		Long	301785	301911	301797	301923	301809	307631	301821	301935	301833	301845	313638
	PLUS	Normal	301874	-	301884	-	307581	-	-	-	-	301894	-
		Long	301879	-	301889	-	307586	-	-	-	-	301899	-

●●● Closed toe ○ Open toe

Product codes

Tights (panty)



			Class 2										
			Beige One (110)		Beige Two (120)		Beige Three (130)		Beige Four (140)		Beige Five (150)	Black	Midnight Blue
			●	○	●	○	●	○	●	○	●	●	●
S	STD	Normal	309863	312564	309873	312526	309883	312549	309893	312538	309903	309913	313644
		Long	309868	312521	309878	312533	309888	312558	309898	312544	309908	309918	313649
	PLUS	Normal	312282	-	312269	-	312292	-	-	-	-	312304	-
		Long	312287	-	312274	-	312297	-	-	-	-	312309	-
M	STD	Normal	309864	312517	309874	312529	309884	312550	309894	312539	309904	309914	313645
		Long	309869	312522	309879	312534	309889	312559	309899	312545	309909	309919	313650
	PLUS	Normal	312283	-	312270	-	312293	-	-	-	-	312305	-
		Long	312288	-	312275	-	312298	-	-	-	-	312310	-
L	STD	Normal	309865	312518	309875	312530	309885	312551	309895	312540	309905	309915	313646
		Long	309870	312523	309880	312535	309890	312560	309900	312546	309910	309920	313651
	PLUS	Normal	312284	-	312271	-	312294	-	-	-	-	312306	-
		Long	312289	-	312276	-	312299	-	-	-	-	312311	-
XL	STD	Normal	309866	312519	309876	312531	309886	312552	309896	312542	309906	309916	313647
		Long	309871	312524	309881	312536	309891	312561	309901	312547	309911	309921	313652
	PLUS	Normal	312285	-	312272	-	312295	-	-	-	-	312307	-
		Long	312290	-	312277	-	312300	-	-	-	-	312312	-
XXL	STD	Normal	309867	312520	309877	312532	309887	312553	309897	312543	309907	309917	313648
		Long	309872	312525	309882	312537	309892	312562	309902	312548	309912	309922	313653
	PLUS	Normal	312286	-	312273	-	312296	-	-	-	-	312308	-
		Long	312291	-	312278	-	312301	-	-	-	-	312313	-

●●● Closed toe ○ Open toe



Mosaic pattern




Chequered pattern



Thigh in Polka dot



Below knee in Chequered



Tights (panty) in Mosaic

Material composition

Chequered and Mosaic

74% Polyamide, 25% Elastane, 1% Cotton (content can vary slightly by product style)

Polka Dot

40% Polyamide Aquarius, 32% Polyamide, 27% Elastane, 1% Polyester

All garments latex free



Scan here for this product's measuring form

Style PATTERNS

A stylish, semi-transparent compression garment designed with breathable, moisture-wicking fabric, soft materials and an easy glide-on design for all-day comfort and compliance.

Compression class

15-20 mmHg

AFNOR 2

Product benefits



Patterns



Semi-transparent



Stylish

Suitable for

- Venous insufficiency
- Mild lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Lipoedema (lipalgia)
- Mixed disease

EMIS search tips

Product name

Style PATTERNS

SIG STY PAT
EMIS TIP

Compression class (mmHg level)

Class 2 15
EMIS TIP

Sizing chart

All measurements are shown in cm

Below knee

	S	M	L	XL
Ankle (cB)	19-22	22-24	24-26	26-29
Calf (cC)	32-38	35-41	38-44	41-47
EMIS TIP	SMA	MED	LA	LA


Thigh with grip top / tights (panty)

	S	M	L	XL
Ankle (cB)	19-22	22-24	24-26	26-29
Calf (cC)	32-38	35-41	38-44	41-47
Thigh (cG)	45-58	48-61	51-65	54-68
Hip (cH) (For tights only)	≤ 105	≤ 115	≤ 125	≤ 125
EMIS TIP	SMA	MED	LA	LA


Length

	Normal	Long
Below knee (ℓD)	≤ 40	≥ 40
Thigh (ℓG)	≤ 75	≥ 75
Tights (panty) (ℓG)	≤ 75	≥ 75
Shoe size	36/40 - 3.5/6.5	37/41 - 4/7.5
EMIS TIP	NO	LO


Below knee



Thigh



Tights (panty)



Styles and colours

● Closed toe

		Class 2		
		Chequered	Mosaic	Polka Dot
Below knee	BEL	●	-	●
Thigh	THI	●	-	●
Tights (panty)	TIG	●	●	●
EMIS TIP		CHE	MOS	DOT



Scan here for this product's ordering codes

Style PATTERNS



Product codes


Below knee



		Class 2	
		Chequered	Polka Dot
		●	●
S	Normal	316856	316832
	Long	316860	316836
M	Normal	316857	316833
	Long	316861	316837
L	Normal	316858	316834
	Long	316862	316838
XL	Normal	316859	316835
	Long	316863	316839

● Closed toe


Thigh



		Class 2	
		Chequered	Polka Dot
		●	●
S	Normal	316864	316840
	Long	316868	316844
M	Normal	316865	316841
	Long	316869	316845
L	Normal	316866	316842
	Long	316870	316846
XL	Normal	316867	316843
	Long	316871	316847

● Closed toe

Tights



		Class 2		
		Chequered	Mosaic	Polka Dot
		●	●	●
S	Normal	316872	316880	316848
	Long	316876	316884	316852
M	Normal	316873	316881	316849
	Long	316877	316885	316853
L	Normal	316874	316882	316850
	Long	316878	316886	316854
XL	Normal	316875	316883	316851
	Long	316879	316887	316855

● Closed toe



Open toe
on below knee
and thigh



SensinnoV
grip top band
on thigh



Below knee
in Black

Thigh with
grip top in
Sable

Material composition

Below knee

32% Polyamide, 25% Cotton, 22% Modal, 21% Elastane

Thigh

36% Elastane, 33% Polyamide, 16% Cotton, 15% Modal

All garments latex free



Scan here
for this product's
measuring form

Essential COTON


Offers year-round comfort with a breathable cotton blend that regulates skin temperature, wicks away moisture and is gentle on sensitive skin.

Compression class



AFNOR 2 AFNOR 3

Product benefits



Suitable for

- Venous insufficiency
- Mild lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Lipoedema (lipalgia)
- Mixed disease
- Sensitive skin

EMIS search tips

Product name
Essential COTON

SIG ESS COT
EMIS TIP

Compression class (mmHg level)

Class 2 15 Class 3 20
EMIS TIP EMIS TIP

Sizing chart

All measurements are shown in cm


Below knee and thigh with grip top

	S	M	L	XL	XXL
Ankle (cB)	19-22	22-24	24-26	26-29	29-32
Calf (cC)	32-39	35-42	38-45	41-48	44-51
Thigh (cG)	46-58	50-62	54-66	58-70	62-74
EMIS TIP	SMA	MED	LAR	LAR	LAR

Length

	Normal	Long
Below knee (lD)	≤ 42	≥ 42
Thigh (lG)	≤ 79	≥ 79
Shoe size	40/43 - 6/9	43/46 - 9/11
EMIS TIP	NOR	LON

Below knee Thigh



Styles and colours

- Closed toe CLO
- Open toe OPE

		Class 2			Class 3	
		Black	Granite	Sable	Black	Sable
Below knee	BEL	● ○	●	● ○	●	●
Thigh with grip top	THI	●	-	● ○	●	●
EMIS TIP		BLA	GRA	SAB	BLA	SAB



Scan here
for this
product's
ordering
codes

Essential COTON



Product codes

Below knee




		Class 2					Class 3	
		Black		Granite	Sable		Black	Sable
		●	○	●	●	○	●	●
S	Normal	303146	303186	303176	303156	303196	303246	303256
	Long	303151	303191	303181	303161	303201	303251	303261
M	Normal	303147	303187	303177	303157	303197	303247	303257
	Long	303152	303192	303182	303162	303202	303252	303262
L	Normal	303148	303188	303178	303158	303198	303248	303258
	Long	303153	303193	303183	303163	303203	303253	303263
XL	Normal	303149	303189	303179	303159	303199	303249	303259
	Long	303154	303194	303184	303164	303204	303254	303264
XXL	Normal	303150	303190	303180	303160	303200	303250	303260
	Long	303155	303195	303185	303165	303205	303255	303265

●●● Closed toe ○● Open toe

Product codes

Thigh with grip top



		Class 2			Class 3	
		Black	Sable		Black	Sable
		●	●	○	●	●
S	Normal	303206	303216	303236	303266	303276
	Long	303211	303221	303241	303271	303281
M	Normal	303207	303217	303237	303267	303277
	Long	303212	303222	303242	303272	303282
L	Normal	303208	303218	303238	303268	303278
	Long	303213	303223	303243	303273	303283
XL	Normal	303209	303219	303239	303269	303279
	Long	303214	303224	303244	303274	303284
XXL	Normal	303210	303220	-	303270	303280
	Long	303215	303225	-	303275	303285

●● Closed toe ○ Open toe



Sensinnov
grip top band
on thigh



**Below
knee in
Black**



**Thigh with grip
top in Black**

Material composition

Below knee

84% Polyamide, 16%
Elastane

Thigh

82% Polyamide, 18%
Elastane

All garments latex free



Scan here
for this product's
measuring form

Essential MICROFIBRE

Delivers exceptional softness with 100% microfibre next to the skin, offering outstanding breathability, moisture control and a secure fit with the Sensinnov grip top band on thigh styles.

Compression class

15-20
mmHg

AFNOR 2

Product benefits



Microfibre Resistance Stylish

Suitable for

- Venous insufficiency
- Mild lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Lipoedema (lipalgia)
- Mixed disease

EMIS search tips

Product name

Essential MICROFIBRE

SIG ESS MIC

EMIS TIP

Sizing chart

All measurements are shown in cm

Below knee

	XS		S			M			L			XL			XXL	
Ankle (cB)	17-19		19-22			22-24			24-26			26-29			29-32	
	SLIM	STD	SLIM	STD	PLUS	SLIM	STD	PLUS	SLIM	STD	PLUS	SLIM	STD	PLUS	STD	PLUS
Calf (cC)	21-28	28-35	24-31	31-38	38-45	27-34	34-41	41-48	29-36	36-43	43-50	32-39	39-46	46-53	42-49	49-56
EMIS TIP	SMA	SMA	SMA	SMA	SMA	MED	MED	MED	LAR	LAR	LAR	LAR	LAR	LAR	LAR	LAR

Thigh with grip top


	S		M		L		XL		XXL
Ankle (cB)	19-22		22-24		24-26		26-29		29-32
	SLIM	STD	SLIM	STD	SLIM	STD	SLIM	STD	STD
Calf (cC)	28-35	31-38	31-38	34-41	33-40	36-43	36-43	39-46	42-49
Thigh (cG)	38-48	48-61	42-52	52-65	47-57	57-70	51-61	61-74	63-76
EMIS TIP	SMA	SMA	MED	MED	LAR	LAR	LAR	LAR	LAR

Length

	Normal	Long
Below knee (lD)	≤ 42	≥ 42
Thigh (lG)	≤ 79	≥ 79
Shoe size	40/43 - 6/9	43/46 - 9/11
EMIS TIP	NOR	LON

Below knee

Thigh



Styles and colours

	EMIS TIP	Class 2	
		Black	Marine*
Below knee	BEL	●	●
Thigh with grip top	THI	●	-
	EMIS TIP	BLA	MAR

* Marine colour only available in standard sizes.



Scan here
for this
product's
ordering
codes

Essential MICROFIBRE



Product codes


Below knee



Class 2	
Black	Marine
●	●

XS	SLIM	Normal	303048	-
		Long	303053	-
	STD	Normal	303000	303036
		Long	303006	303042
S	SLIM	Normal	303049	-
		Long	303054	-
	STD	Normal	303001	303037
		Long	303007	303043
	PLUS	Normal	303068	-
		Long	303073	-
M	SLIM	Normal	303050	-
		Long	303055	-
	STD	Normal	303002	303038
		Long	303008	303044
	PLUS	Normal	303069	-
		Long	303074	-
L	SLIM	Normal	303051	-
		Long	303056	-
	STD	Normal	303003	303039
		Long	303009	303045
	PLUS	Normal	303070	-
		Long	303075	-
XL	SLIM	Normal	303052	-
		Long	303057	-
	STD	Normal	303004	303040
		Long	303010	303046
	PLUS	Normal	303071	-
		Long	303076	-
XXL	STD	Normal	303005	303041
		Long	303011	303047
	PLUS	Normal	303072	-
		Long	303077	-

Thigh with grip top



Class 2	
Black	
●	

S	SLIM	Normal	303120
		Long	303124
	STD	Normal	303100
		Long	303105
M	SLIM	Normal	303121
		Long	303125
	STD	Normal	303101
		Long	303106
L	SLIM	Normal	303122
		Long	303126
	STD	Normal	303102
		Long	303107
XL	SLIM	Normal	303123
		Long	303127
	STD	Normal	303103
		Long	303108
XXL	STD	Normal	303104
		Long	303109

● Closed toe

●● Closed toe



Below knee in Black / Aubergine

Below knee in Blue / Grey

Material composition

Below knee

38% Polyamide, 36% Cotton, 15% Polypropylene, 11% Elastane
Latex free



Scan here for this product's measuring form

Style COLOURS

Offers a soft feel against the skin, with a low-profile toe seam and an ergonomically designed sole that cushions the ball of the foot for all-day comfort.

Compression class

15-20 mmHg

AFNOR 2

Product benefits



Moisture-wicking Stylish Cushioned sole

Suitable for

- Venous insufficiency
- Mild lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Mixed disease

EMIS search tips

Product name

Style COLOURS

SIG STY COL

EMIS TIP

Sizing chart

All measurements are shown in cm


Below knee

	S	M	L	XL	XXL
Ankle (cB)	19-22	22-24	24-26	26-29	29-32
Calf (cC)	31-38	34-41	36-43	39-46	42-49
EMIS TIP	SMA	MED	LAR	LAR	LAR

Length

	Normal	Long
Below knee (lD)	≤ 42	≥ 42
Shoe size	40/43 - 6/9	43/46 - 9/11
EMIS TIP	NOR	LON

Below knee



Styles and colours

● Closed toe

		Class 2	
		Black / Aubergine	Blue / Grey
Below knee	EMIS TIP	BEL	BEL
	EMIS TIP	BLA	BLU



Scan here for this product's ordering codes

Style COLOURS



Product codes

Below knee



		S		M		L		XL		XXL	
		Normal	Long	Normal	Long	Normal	Long	Normal	Long	Normal	Long
Black / Aubergine	●	314943	314948	314944	314949	314945	314950	314946	314951	314947	314952
Blue / Grey	●	314953	314958	314954	314959	314955	314960	314956	314961	314957	314962



Below knee in Carbon



Below knee in Black



Below knee in Cyanide

Material composition

Below knee

81% Polyamide, 19% Elastane
Latex free



Scan here for this product's measuring form

Active MASCULINE

Features discreet ribbed styling, hydroregulating fibres, and built-in zinc odour protection, with cushioned padding and anatomical support for all-day comfort.

Compression class

23-32 mmHg

RAL 2

Product benefits



Moisture-wicking



Stylish



Cushioned sole

Suitable for

- Venous insufficiency
- Mild Lymphoedema / chronic oedema
- Prevention of ulcer recurrence
- Lipoedema (lipalgia)
- Deep vein thrombosis (DVT)

EMIS search tips

Product name

Active MASCULINE

SIG ACT MAS

EMIS TIP

Sizing chart

All measurements are shown in cm


Below knee

	XXS		XS		S		M		L		XL	
Ankle (cB)	18-20		20-22		22-24		24-26		26-29		29-32	
	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS
Calf (cC)	25-32	29-36	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49	42-49	43-51
EMIS TIP	SMA	SMA	SMA	SMA	MED	MED	LAR	LAR	LAR	LAR	LAR	LAR

Length

	Normal	Long
Below knee (cD)	35-40	41-45
EMIS TIP	NOR	LON

Below knee



Styles and colours

●●● Closed toe

		Class 2		
		Black	Carbon	Cyanide
Below knee	EMIS TIP	●	●	●
	EMIS TIP	BLA	GRA	CYA



Scan here for this product's ordering codes



Product codes

Below knee



			Class 2		
			Black	Carbon	Cyanide
			●	●	●
XXS	STD	Normal	72060	72012	72036
		Long	72061	72013	72037
	PLUS	Normal	72072	72024	72048
		Long	72073	72025	72049
XS	STD	Normal	72062	72014	72038
		Long	72063	72015	72039
	PLUS	Normal	72074	72026	72050
		Long	72075	72027	72051
S	STD	Normal	72064	72016	72040
		Long	72065	72017	72041
	PLUS	Normal	72076	72028	72052
		Long	72077	72029	72053
M	STD	Normal	72066	72018	72042
		Long	72067	72019	72043
	PLUS	Normal	72078	72030	72054
		Long	72079	72031	72055
L	STD	Normal	72068	72020	72044
		Long	72069	72021	72045
	PLUS	Normal	72080	72032	72056
		Long	72081	72033	72057
XL	STD	Normal	72070	72022	72046
		Long	72071	72023	72047
	PLUS	Normal	72082	72034	72058
		Long	72083	72035	72059

●●● Closed toe




Made to measure (MTM)

Products available

MTM Essential COMFORTABLE	48-49
MTM Essential THERMOREGULATING	50-51
MTM TRADITIONAL	52-53
MTM Style SEMITRANSSPARENT	54-55
MTM OPTIFORM Flex - Flat knit	56-57
MTM OPTIFORM Hold - Flat knit	



For all Made to Measure (MTM) garments measuring points, please use the **Circular or Flat Knit Made to Measure form** or our **Made to Measure brochure** (available via the QR code).



Knobbed grip
top band



Sensinnov
grip top band
on thigh

Material composition

All styles

65% Polyamide, 35%
Elastane (content
can vary slightly by
compression choice)

Latex free



Scan here
for the MTM Circular
knit measuring form

MTM Essential COMFORTABLE

The modern fabric provides a soft and smooth feel against the skin while offering breathability and moisture-wicking to aid skin care.

Compression class

18-21
mmHg

23-32
mmHg

RAL 1

RAL 2

Product benefits



Soft



Easy
donning &
doffing



Mois-
ture-
wicking

Suitable for

- Venous insufficiency
- Lymphoedema / chronic oedema
- Sensitive skin
- Regulating skin temperature
- Active leg ulcers
- Mixed disease
- Lipoedema (lipalgia)
- Prevention of leg ulcers and oedema



Below knee in Black
with Closed toe



**Thigh with waist
attachment** in Savannah
with Open toe



**Thigh with
grip top** in
Grey with
Closed toe

Tights (panty) in Skin
with Closed toe

EMIS search tips

Product name

MTM Essential COMFORTABLE

MTM ESS COMF
EMIS TIP

Compression class (mmHg level)

Class 1 18
EMIS TIP

Class 2 23
EMIS TIP

Styles and colours

Below knee






Thigh

Thigh w/attachment

Tights (panty)

Tights (panty)
one leg

Tights (panty)
maternity



Left and right

Left and right

●●●● Closed toe
○○○○ Open toe

CLO
OPE
EMIS TIP

		Black	Caramel	Grey	Savannah	Skin
Below knee	BEL	● ○	● ●	● ○	● ●	● ●
Below knee with grip top	BEL GRI	● ○	● ●	● ○	● ●	● ●
Thigh without grip top	THI	● ○	● ●	● ○	● ●	● ●
Thigh with grip top	THI GRI	● ○	● ●	● ○	● ●	● ●
Thigh waist attachment	WAI LEF	● ○	● ●	● ○	● ●	● ●
	WAI RIG	● ○	● ●	● ○	● ●	● ●
Tights (panty)	PAN	● ○	● ●	● ○	● ●	● ●
Tights (panty) bodyform	BOD	● ○	● ●	● ○	● ●	● ●
Tights (panty) one leg	ONE	● ○	● ●	● ○	● ●	● ●
Tights (panty) maternity*		● ○	● ●	● ○	● ●	● ●

* Product not available on FP10 / prescription.



Scan here
to view our
MTM Circular
knit brochure



Knobbed grip top band



Sensinno grip top band on thigh

Material composition

All styles

52% Polyamide, 34% Elastane, 14% cotton (content can vary slightly by compression choice)

Latex free



Scan here for the MTM Circular knit measuring form

MTM Essential THERMOREGULATING

Combines all day comfort in a classic design with 100% cotton against the skin. The thermoregulating properties make it ideal for any season.

Compression class

18-21 mmHg RAL 1 23-32 mmHg RAL 2 34-46 mmHg RAL 3

Product benefits

Sensitive skin 100% cotton against the skin Breathable

Suitable for

- Venous insufficiency
- Sensitive skin
- Regulating skin temperature
- Active leg ulcers
- Mixed disease
- Lipoedema (lipalgia)
- Prevention of leg ulcers and oedema



Below knee in Black with Open toe



Thigh with grip top in Black with Open toe



Thigh with waist attachment in Nature with Closed toe



Tights (panty) in Nature with Closed toe

EMIS search tips

Product name

MTM Essential THERMOREGULATING

MTM ESS THERMO EMIS TIP

Compression class (mmHg level)

Class 1 18 EMIS TIP Class 2 23 EMIS TIP Class 3 34 EMIS TIP

Styles and colours

Below knee Thigh Thigh w/attachment Tights (panty) Tights (panty) one leg Tights (panty) maternity

Left and right Left and right

	EMIS TIP	Toe	
		Black	Nature
Below knee	BEL	● ○	● ○
Below knee with grip top	BEL GRI	● ○	● ○
Thigh without grip top	THI	● ○	● ○
Thigh with grip top	THI GRI	● ○	● ○
Thigh waist attachment	WAI LEF	● ○	● ○
	WAI RIG	● ○	● ○
Tights (panty)	PAN	● ○	● ○
Tights (panty) bodyform	BOD	● ○	● ○
Tights (panty) one leg	ONE	● ○	● ○
Tights (panty) maternity*		● ○	● ○

* Product not available on FP10 / prescription.



Scan here to view our MTM Circular knit brochure



Knobbed grip top band

Material composition

All styles

70% Polyamide, 30% Natural rubber latex (content can vary slightly by compression choice)

Latex free



Scan here for the MTM Circular knit measuring form

MTM TRADITIONAL 500

Offers highly effective compression that is durable and comfortable for those with more complex indications.

Compression class

23-32 mmHg RAL 2 34-46 mmHg RAL 3 >49 mmHg RAL 4

Product benefits

Natural rubber latex Durable Containment

Suitable for

- Severe venous conditions
- Chronic oedema
- Lymphoedema
- Prevention of ulcer recurrence



Thigh with waist attachment



Tights (panty)



Thigh with grip top

EMIS search tips

Product name

MTM TRADITIONAL 500

MTM TRAD
EMIS TIP

Compression class (mmHg level)

Class 2 23
EMIS TIP

Class 3 34
EMIS TIP

Class 4 49
EMIS TIP

Styles and colours

Below knee Half Thigh Thigh Thigh w/ attachment Tights (panty) Tights (panty) one leg Tights (panty) maternity

Left and right Left and right

● Closed toe CLO
● Open toe OPE
EMIS TIP


Beige

Below knee	BEL	● ●	
Below knee with grip top	BEL GRI	● ●	
Half thigh*		● ●	
Half thigh with grip top*		● ●	
Thigh without grip top	THI	● ●	
Thigh with grip top	THI GRI	● ●	
Thigh waist attachment	WAI LEF	● ●	
	WAI RIG	● ●	
Tights (panty)	PAN	● ●	
Tights (panty) one leg	PAN	● ●	
Tights (panty) bodyform	BOD	● ●	
Tights (panty) maternity*		● ●	
Full compression	Half thigh*	● ●	
	Half thigh with grip top*	● ●	
	Thigh without grip top	THI FULL	● ●
	Thigh with grip top	THI GRI FULL	● ●
	Thigh waist attachment	WAI LEF FULL	● ●
		WAI RIG FULL	● ●
	Tights (panty)	PAN FULL	● ●
	Tights (panty) one leg	PAN FULL	● ●
Tights (panty) maternity*		● ●	

* Product not available on FP10 / prescription.



Scan here to view our MTM Circular knit brochure



Knobbed
grip top band on
thigh



Sensinno
laced grip top
on thigh



*Thigh with
grip top in
Marine*



*Below knee
in Grey with Open toe*



*Tights (panty) in
Savannah*

Material composition

Below knee

61% Polyamide, 39% Elastane

Thigh

63% Polyamide,
37% Elastane

Tights

59% Polyamide,
41% Elastane

All garments latex free



Scan here
for the MTM Circular
knit measuring form

MTM Style SEMITRSPARENT

Stylish and in a range of colours, Style SEMITRSPARENT fits perfectly to the shape of the leg for an elegant look.

Compression class

18-21
mmHg

RAL 1

23-32
mmHg

RAL 2

Product benefits



Semi-transparent



Durable



Breathable

Suitable for

- Venous insufficiency
- Mixed disease
- Lipalgia (previously lipoedema)
- Prevention of leg ulcers and oedema
- Pre & post vein treatment with minimal oedema

EMIS search tips

Product name

MTM Style SEMITRSPARENT

MTM STY SEMI
EMIS TIP


Compression class (mmHg level)

Class 1 18
EMIS TIP


Class 2 23
EMIS TIP

Styles and colours


Below knee




Thigh



Tights (panty)



Tights (panty) maternity



●●●●● Closed toe

CLO

○●●●○ Open toe

OPE

EMIS TIP

		Black	Caramel	Grey	Marine	Savannah	Skin
Below knee	BEL	● ○	● ○	● ○	● ○	● ○	● ○
Below knee with grip top	BEL GRI	● ○	● ○	● ○	● ○	● ○	● ○
Thigh with grip top	THI GRI	● ○	● ○	● ○	● ○	● ○	● ○
Tights (panty)	PAN	● ○	● ○	● ○	● ○	● ○	● ○
Tights (panty) maternity*		● ○	● ○	● ○	● ○	● ○	● ○

* Product not available on FP10 / prescription.



Scan here
to view our
MTM Circular
knit brochure

Which to choose?

FLEX

Strong containment but more stretch and freedom of movement

HOLD

Stronger containment with higher static stiffness

Material composition

FLEX - All styles

30% Polyamide, 25% Elastane, 20% Cotton, 19% Modal, 6% Viscose (content can vary slightly by compression choice)

HOLD - All styles

64% Polyamide, 20% Cotton, 16% Elastane (content can vary slightly by compression choice)

Both materials are latex free



Scan here for the MTM OPTIFORM flatknit measuring form

MTM OPTIFORM Flat knit

Offers highly effective compression that is durable and comfortable for those with more complex indications.

Compression class

18-21 mmHg 23-32 mmHg 34-46 mmHg

RAL 1 RAL 2 RAL 3*

*RAL 3 only available in OPTIFORM HOLD

Product benefits



Reinforced zones Containment Made to measure

Suitable for

- Venous Insufficiency
- Lymphoedema/ chronic oedema
- Lipoedema (lipalgia)
- Mixed Disease
- Leg Ulcer prevention

Tights (panty)




Straight toe
Available for Closed and Open toe



Slant toe
Available for Closed and Open toe

Styles and colours




- Closed toe
- Open toe
- ⊗ No foot

	OPTIFORM Flex					
	OPTIFORM Hold					
	Class 1		Class 2		Class 3	
	Beige	Black	Beige	Black	Beige	Black
Ankle sock	● ○	● ○	● ○	● ○	● ○	● ○
Below knee	● ○	● ○	● ○	● ○	● ○	● ○
Below knee sleeve	⊗	⊗	⊗	⊗	⊗	⊗
Thigh	● ○	● ○	● ○	● ○	● ○	● ○
Thigh sleeve	⊗	⊗	⊗	⊗	⊗	⊗
Thigh with waist attachment (left and right)	● ○	● ○	● ○	● ○	● ○	● ○
Body bandage	⊗	⊗	⊗	⊗	⊗	⊗
Tights (panty)	● ○	● ○	● ○	● ○	● ○	● ○
Tights (panty) with one leg	● ○	● ○	● ○	● ○	● ○	● ○
Shorts	⊗	⊗	⊗	⊗	⊗	⊗
Capri length	⊗	⊗	⊗	⊗	⊗	⊗
Leggings	⊗	⊗	⊗	⊗	⊗	⊗

Accessories



- Silicone strip**
Silicone strips can be added anywhere on the garment to help keep it in place during wear.
- Knobbed grip top**
Knobbed grip top placed at the top of the garment. Helps to hold the garment in place during wear.
- Knee functional zone**
Enhances maximum freedom of movement of the joint.
- Adjustable waistband**
A light band inside the waistband, which is adjustable and allows the wearer to adjust the band throughout the day for additional comfort.
- Liner pocket**
Can be used to protect the skin from friction especially at the popliteal area or can be used to hold foam pads to treat troublesome fibrotic tissues on the limb.
- Ankle pad (1 piece)**
These are applied to the garment at the malleolar area of the ankle. Provides extra pressure / support and help to shape the area, reducing swelling.



- Fly opening**
This can be used for ease of use and in some cases when pads need to be used with hosiery.
- Reinforced gusset**
Provides extra support at the gusset area when managing genital oedema or for extra support to reduce the risk of genital oedema when wearing compression.
- Zip with lining**
Indicated for use with patients who are unable to pull garments on by themselves. Zip fastening can aid garment application.
- Quad plus**
Designed for those patients who may have larger thigh areas where standard garments may not fit. This allows more fabric to fit and support the leg shape.



Scan here to view our MTM OPTIFORM Flatknit brochure



Inelastic compression wraps

Products available

COMPREFLEX Standard - Calf	62-65
COMPREFLEX Accessories	66
COOLFLEX Standard - Calf	68-69
COMPREFLEX Standard - Thigh	70-71
COMPREFLEX Standard - Knee	72-73
COMPREBOOT Standard	74-75

It's a wrap



Burnout

High nurse-patient workloads are the main cause of burnout among NHS nurses.² 61% of nurses felt too busy to provide the level of care they desired.³



Caseload pressure

District nurse caseloads over 600 patients rose from 11.5% in 2019 to 16.2% in 2023. More admin, travel and virtual ward responsibilities compound the load.⁷



Reduce application time

Using inelastic compression wraps takes typically 5-10 minutes to apply, a 66% reduction compared to traditional 4-layer bandaging.^{4,5}



Reduce visits

Treating patients with inelastic compression wraps gives patients the ability to self-care, with the potential to reduce nurse visits.⁶

Empowering Collaborative Care

Sigvaris inelastic compression wraps offer 3 adjustable compression levels from 20 to 50 mmHg of compression in 1 wrap, giving the ability to measure pressure and adjust treatment.

COMPREFLEX Calf

AccuTabs® – Patented adjustable pressure monitoring tabs for easy, accurate and consistent control

Low profile – Designed to fit discreetly under clothing, including trousers

Fast drying – Machine washable at 30°, dries quickly*

Premium materials – Soft, breathable Breathe-O-Prene® offers a comfortable fit. Moisture-regulating and gentle on the skin



COOLFLEX Standard - Calf

Precise pressure adjustment – Patented easy-to-use hook and eye system for accurate pressure adjustment

Stay cool – Lightweight, breathable mesh keeps skin comfortable

Low-profile – Slim design fits easily under clothing and with footwear

Fast drying – Machine washable at 30°, dries quickly*

*Drying time may vary depending on conditions.



COMPREFLEX Transition Calf
in Black

COMPREFLEX Complete
in Beige

COMPREFLEX Standard - Calf
in Black

AccuTab® allows the patient and HCP to accurately measure the compression required

Material composition

37% Polyamide, 27% Polyurethane Foam, 26% Cotton, 6% Elastane, 4% Polyester
Latex free



Scan here for this product's measuring form

COMPREFLEX Calf

Featuring a contoured calf for a superior fit, the wrap offers adjustable compression with AccuTabs®, easy application with a front stretch panel, and durable double-stitched composition.

Compression class

20-50 mmHg

Variable

Product benefits



Containment Durable

Easy donning & doffing

Suitable for

- Lymphoedema
- Oedema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including ulcer due to venous stasis
- Lipoedema

EMIS search tips

Product name
COMPREFLEX Calf

SIG COMPREFL
EMIS TIP

Compression class (mmHg level)
Variable 20
EMIS TIP

Sizing chart

All measurements are shown in cm

Calf


	S	M	L	XL	XXL
Ankle (B)	19-26	21-30	26-36	31-41	36-46
Calf (C)	29-39	34-44	40-50	46-56	54-64
EMIS TIP	SMA	MED	LAR	LAR	LAR

Length

	Regular	Tall
Below knee (L1)	30-33.9	34-42
EMIS TIP	SHO	LON

For leg length measurements, take from outer ankle bone to fibular head.

Below knee



Measure the points at the heights given
E.G.: Circumference for B is measured 5cm above the outer ankle bone

Styles and colours

		Beige	Black
COMPREFLEX Standard - Calf	BEL	●	●
COMPREFLEX Transition - Calf	CAL TRANS	●	●
COMPREFLEX Complete - Calf	CAL COMPL	●	●
EMIS TIP		BEI	BLA



Scan here for this product's ordering codes

COMPREFLEX
Standard Calf

20-50
mmHg

Variable

Features and benefits

- Non-compressive liner
- Soft, conforming and breathable fabric provides a comfortable fit and wicks away moisture
- Low profile under clothing
- Available in Black or Beige
- Strap extenders may be applied if the patient exceeds standard size ranges



Included with this product



AccuTabs® allows the patient and HCP to accurately measure the compression required



Basic Liner is non compressive

Basic liner

Features and benefits

- Non-compressive liner
- Easy to apply knee-high liners
- *drilease®* moisture-wicking and odour control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties

Product materials

- 74% Cotton, 13% Nylon, 10% Polyester, 3% Elastane

COMPREFLEX
Transition Calf

20-50
mmHg

Variable

Features and benefits

- Transition Liner provides 10-15mmHg in the foot and ankle
- Soft, conforming and breathable fabric provides a comfortable fit and wicks away moisture
- Easily adjustable
- Available in Black or Beige
- Strap extenders may be applied if the patient exceeds standard size ranges



Included with this product



AccuTabs® allows the patient and HCP to accurately measure the compression required



Transition Liner gives 10-15 mmHg of mild compression in the foot and ankle

Transition liner

Features and benefits

- Mild compression in foot and ankle only
- Comfortable knee-high liners
- *drilease®* moisture-wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

Product materials

- 67% Cotton, 19% Polyamide, 10% Polyester, 4% Elastane

COMPREFLEX
Complete Calf

20-50
mmHg

Variable

Features and benefits

- The Complete Liner provides 10-15mmHg in the calf, foot and ankle
- Soft, conforming and breathable fabric provides a comfortable fit and wicks away moisture
- Easily adjustable
- Available in Black or Beige
- Strap extenders may be applied if the patient exceeds standard size ranges



Included with this product



AccuTabs® allows the patient and HCP to accurately measure the compression required



Complete Liner gives 10-15 mmHg of mild compression in the calf, foot and ankle

Complete liner

Features and benefits

- Mild compression in calf, foot and ankle with graduated compression in calf
- Comfortable knee-high liners
- *drilease®* moisture-wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

Product materials

- 66% Cotton, 19% Polyamide, 9% Polyester, 6% Elastane

COMPREFLEX Calf



Product codes

Below knee



COMPREFLEX Standard Calf

		S	M	L	XL	XXL
Beige	Regular	1411-NF-BKR	1412-NF-BKR	1413-NF-BKR	1414-NF-BKR	1415-NF-BKR
	Tall	1411-NF-BKT	1412-NF-BKT	1413-NF-BKT	1414-NF-BKT	1415-NF-BKT
Black	Regular	1401-NF-BKR	1402-NF-BKR	1403-NF-BKR	1404-NF-BKR	1405-NF-BKR
	Tall	1401-NF-BKT	1402-NF-BKT	1403-NF-BKT	1404-NF-BKT	1405-NF-BKT



COMPREFLEX Transition Calf

		S	M	L	XL	XXL
Beige	Regular	1411-UC-BKR	1412-UC-BKR	1413-UC-BKR	1414-UC-BKR	1415-UC-BKR
	Tall	1411-UC-BKT	1412-UC-BKT	1413-UC-BKT	1414-UC-BKT	1415-UC-BKT
Black	Regular	1401-UC-BKR	1402-UC-BKR	1403-UC-BKR	1404-UC-BKR	1405-UC-BKR
	Tall	1401-UC-BKT	1402-UC-BKT	1403-UC-BKT	1404-UC-BKT	1405-UC-BKT



COMPREFLEX Complete Calf

		S	M	L	XL	XXL
Beige	Regular	1411-FC-BKR	1412-FC-BKR	1413-FC-BKR	1414-FC-BKR	1415-FC-BKR
	Tall	1411-FC-BKT	1412-FC-BKT	1413-FC-BKT	1414-FC-BKT	1415-FC-BKT
Black	Regular	1401-FC-BKR	1402-FC-BKR	1403-FC-BKR	1404-FC-BKR	1405-FC-BKR
	Tall	1401-FC-BKT	1402-FC-BKT	1403-FC-BKT	1404-FC-BKT	1405-FC-BKT



Cotton liner

Features and benefits

- Non-compressive liner
- Made from cotton, it is gentle on the skin and breathable, making it ideal for individuals with sensitive skin.
- Available for calf (one size) and thigh (S / M / L / XL)



EMIS TIP

SIG COTT
LIN WHI


AccuTab® Device*



4 pack

Easy to use

- Position the COMPREFLEX device on the limb, gently position the straps without stretching or applying compression
- Attach one AccuTab® at the end of each strap
- Pull each strap to the prescribed pressure indicated by the lines on the AccuTab®: 20–30, 30–40, or 40–50 mmHg



* AccuTabs® are provided with the COMPREFLEX garment

Strap extenders

Features and benefits

- Provides up to 10cm of additional circumferential size range to COMPREFLEX wraps
- Fixes securely
- Potential solution for patients outside of standard size ranges of COMPREFLEX wraps

Styles and colours

	EMIS TIP	Beige	Black
COMPREFLEX Strap extenders	SIG COMPREF STR	●	●
	EMIS TIP	BEI	BLA



COOLFLEX Standard - Calf

- **Breathable mesh textile** made from mesh and biocompatible materials.
- **Lightweight inelastic wrap** designed to be both light and breathable.
- **Moisture regulating** mesh textile allows wearers to feel cool and comfortable.

**Breathable,
light & comfortable**



COMPREFLEX Accessories

Product codes

Basic liner



One size fits all

1700-NK1

Transition liner



Beige	1611-CK	1612-CK	1613-CK	1614-CK	1615-CK
Black	1601-CK	1602-CK	1603-CK	1604-CK	1605-CK

Complete liner



Beige	1611-CK-F1	1612-CK-F1	1613-CK-F1	1614-CK-F1	1615-CK-F1
Black	1601-CK-F1	1602-CK-F1	1603-CK-F1	1604-CK-F1	1605-CK-F1

Cotton liner

	One size fits all
Thigh Length S / M / L	1793-TH
Thigh Length XL	1794-TH
Calf Length One size fits all	1795-BK1

Strap extenders

Beige	1411-SE1
Black	1401-SE1



Material composition

47% Polyester, 25% Polyamide, 20% Cotton, 8% Elastane, 1% Stainless steel

Latex free



Scan here for this product's measuring form

COOLFLEX Standard - Calf

Combining innovation and practicality, the COOLFLEX Standard - Calf's mesh fabric enhances heat and moisture transfer, while stainless steel hooks ensure durability and ease of use.

Compression class

20-50 mmHg

Variable

Product benefits



Containment



Breathable



Moisture-wicking

Suitable for

- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including ulcer due to venous stasis
- Lymphoedema
- Oedema
- Lipoedema

EZ - Close pressure system with visual indicators. Each closure tab features three white markings indicating the compression ranges from 20-50 mmHg.

EMIS search tips

Product name

COOLFLEX Standard - Calf

COOL
EMIS TIP

Compression class (mmHg level)

Variable 20
EMIS TIP

Sizing chart

All measurements are shown in cm

Calf

	S	S Max	M/L	M/L Max	XL/XXL
Ankle (B)	16-26		23-33		30-42
Below Calf (B1)	22-33	27-38	29-43	35-49	37-54
Calf (C)	28-39	35-50	34-50	45-64	45-64
Above Calf (D)	28-39	35-50	34-50	45-64	45-64
EMIS TIP	SMA	SMA MAX	MED	MED MAX	LAR

Length


	Regular	Tall
Below knee (L1)	30-35.9	36-42
EMIS TIP	REG	TAL

For leg length measurements, take from outer ankle bone to fibular head.

Leg

	Left	Right
EMIS TIP	LEF	RIG

Below knee



Scan here for this product's ordering codes

COOLFLEX Standard - Calf




Product codes

Below knee



			S	S Max	M/L	M/L Max	XL/XXL
Black	Regular	Left	2101-NF-BKR-L	2101-NF-BMR-L	2103-NF-BKR-L	2103-NF-BMR-L	2105-NF-BKR-L
		Right	2101-NF-BKR-R	2101-NF-BMR-R	2103-NF-BKR-R	2103-NF-BMR-R	2105-NF-BKR-R
	Tall	Left	2101-NF-BKT-L	2101-NF-BMT-L	2103-NF-BKT-L	2103-NF-BMT-L	2105-NF-BKT-L
		Right	2101-NF-BKT-R	2101-NF-BMT-R	2103-NF-BKT-R	2103-NF-BMT-R	2105-NF-BKT-R



AccuTab® allows the patient and HCP to accurately measure the compression required

COMPREFLEX Standard - Thigh in Black



COMPREFLEX Standard - Thigh in Beige

Material composition

45% Cotton, 28% Polyamide, 14% Polyurethane Foam, 9% Neoprene, 4% Elastane
Latex free



Scan here for this product's measuring form

COMPREFLEX Standard - Thigh

An adjustable inelastic wrap for lymphoedema, featuring AccuTabs®, breathable fabric, and a secure hip attachment for a comfortable, contoured fit.

Compression class

20-50 mmHg

Variable

Product benefits



Containment



Durable



Easy donning & doffing

Suitable for

- Lymphoedema
- Oedema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including ulcer due to venous stasis
- Lipoedema

EMIS search tips

Product name

COMPREFLEX Standard - Thigh

SIG COMPREFL
EMIS TIP

Compression class (mmHg level)

Variable 20
EMIS TIP

Sizing chart


All measurements are shown in cm

Thigh	S	M	L	XL
Above knee (E)	38-48	46-56	53-63	63-73
Thigh (F)	43-53	51-61	58-68	68-78
Below gluteal fold (G)	48-58	56-66	64-74	74-84
EMIS TIP	SMA	MED	LAR	LAR

Length

	Regular	Tall
Above knee (L2)	20-24.9	25-30
EMIS TIP	SHO	LON

Thigh



Styles and colours

		Beige	Black
COMPREFLEX Standard - Thigh Left	THI LEF	●	●
COMPREFLEX Standard - Thigh Right	THI RIG	●	●
EMIS TIP		BEI	BLA




Scan here for this product's ordering codes

COMPREFLEX Standard - Thigh




Product codes

Thigh



			S	M	L	XL
Beige	Regular	Left	1411-TCR-L	1412-TCR-L	1413-TCR-L	1414-TCR-L
		Right	1411-TCR-R	1412-TCR-R	1413-TCR-R	1414-TCR-R
	Tall	Left	1411-TCT-L	1412-TCT-L	1413-TCT-L	1414-TCT-L
		Right	1411-TCT-R	1412-TCT-R	1413-TCT-R	1414-TCT-R
Black	Regular	Left	1401-TCR-L	1402-TCR-L	1403-TCR-L	1404-TCR-L
		Right	1401-TCR-R	1402-TCR-R	1403-TCR-R	1404-TCR-R
	Tall	Left	1401-TCT-L	1402-TCT-L	1403-TCT-L	1404-TCT-L
		Right	1401-TCT-R	1402-TCT-R	1403-TCT-R	1404-TCT-R



AccuTab® allows the patient and HCP to accurately measure the compression required

COMPREFLEX Standard - Knee in Black



COMPREFLEX Standard - Knee in Beige



Material composition

50% Polyamide 42% Polyurethane Foam, 8% Elastane
Latex free



Scan here for this product's measuring form

COMPREFLEX Standard - Knee

With a built-in stretch panel for easy application and a four-way strap for targeted support, this product ensures secure, effective compression in hard-to-reach areas like the patella or popliteal space.

Compression class

20-50 mmHg

Variable

Product benefits



Containment Durable Easy donning & doffing

Suitable for

- Lymphoedema
- Oedema
- Chronic venous disease (CVD)
- Chronic venous insufficiency (CVI) including ulcer due to venous stasis
- Lipoedema

EMIS search tips

Product name

COMPREFLEX Standard - Knee

SIG COMPRE

EMIS TIP


Sizing chart

All measurements are shown in cm

Knee

	S	M	L	XL
Below knee (D)	29-39	34-44	39-49	44-65
Above knee (E)	38-48	46-56	53-63	58-73
EMIS TIP	SMA	MED	LAR	LAR

Knee



Measure the points at the heights given
E.G.: Circumference for E is measured 5cm above top of the patella

Styles and colours

		Beige	Black
COMPREFLEX Standard - Knee	KNE	●	●
	EMIS TIP	BEI	BLA




Scan here for this product's ordering codes

COMPREFLEX Standard - Knee



Product codes

Knee



	S	M	L	XL
Beige	1411-KP	1412-KP	1413-KP	1414-KP
Black	1401-KP	1402-KP	1403-KP	1404-KP



COMPREBOOT Standard in Black



COMPREBOOT Standard and COMPREFLEX Standard - Calf in Beige

Material composition

67% Neoprene,
30% Polyamide, 3%
Elastane
Latex free



Scan here for this product's measuring form

COMPREBOOT Standard

Designed for optimal comfort and effectiveness while allowing you to apply compression in the foot. Its low-stretch material enhances fit and simplifies application, ensuring reliable support.

Compression class

30-40 mmHg

Variable

Product benefits



Containment



Durable



Easy donning & doffing

Suitable for

- Lymphoedema / chronic oedema
- Foot oedema

EMIS search tips

Product name
COMPREBOOT Standard


SIG COMPREB
EMIS TIP

Sizing chart

All measurements are shown in cm

Foot	S		M / L		XL / XXL	
	Regular	Long	Regular	Long	Regular	Long
Ankle bend and heel (Y)	28-36		≤ 39		≤ 44	
Metatarsal heads (A)	22-26		≤ 30		≤ 33	
Inside of foot (K)	14-18	19-23	18-20	23-25	20-22	25-27
EMIS TIP	SMA REG	SMA LON	MED REG	MED LON	LAR REG	LAR LON

Boot



Styles and colours

COMPREBOOT Standard	EMIS TIP	Beige	Black
		BOO	BOO
	EMIS TIP	BEI	BLA



Scan here for this product's ordering codes

COMPRESSBOOT Standard



Product codes

Foot



	S		M / L		XL / XXL	
	Regular	Long	Regular	Long	Regular	Long
Beige	1811-BTR	1811-BTL	1813-BTR	1813-BTL	1815-BTR	1815-BTL
Black	1801-BTR	1801-BTL	1803-BTR	1803-BTL	1805-BTR	1805-BTL



Upper limb

Products available

COMPREFLEX - Standard Arm	78-79
ADVANCE Armsleeve	80-81
TRADITIONAL Armsleeve	82-83



Scan the QR code to view our Upper limb brochure which includes measuring points

COMPREFLEX Standard - Arm in Beige



COMPREFLEX Standard - Arm in Black

Material composition

43% Polyamide, 26% Polyurethane Foam, 20% Cotton, 6% Elastane, 5% Neoprene
Latex free



Scan here for this product's measuring form

COMPREFLEX Standard - Arm

Offers gentle, adjustable compression with a one-handed donning design, breathable comfort and a 90° elbow bend for natural movement.

Compression class

20-30 mmHg

Variable

Product benefits



Containment



Durable



Easy donning & doffing

Suitable for

- Lymphoedema
- Oedema
- Lipoedema
- Night time wear to prevent rebound swelling

EMIS search tips

Product name

COMPREFLEX Standard - Arm

SIG COMPRE
EMIS TIP

Sizing chart

All measurements are shown in cm


Armsleeve

	S	M	L	XL
Wrist (A)	14-19	15-20	16-21	18-23
Elbow (B)	20-27	24-31	28-35	32-39
Mid bicep (B1)	22-30	26-34	30-39	35-44
Axilla (C)	23-32	28-37	33-43	39-49
EMIS TIP	SMA	MED	LAR	LAR

Length

	Short	Regular	Long
Arm (D)	38-42.9	43-47.9	48-53
EMIS TIP	SHO	REG	LON

Armsleeve



Styles and colours

		Beige	Black
COMPREFLEX Standard - Arm	EMIS TIP	ARM	ARM
	EMIS TIP	BEI	BLA



Scan here for this product's ordering codes

COMPREFLEX Standard - Arm



Product codes

Arm



		S	M	L	XL
Beige	Short	1411-ARS	1412-ARS	1413-ARS	1414-ARS
	Regular	1411-ARR	1412-ARR	1413-ARR	1414-ARR
	Long	1411-ARL	1412-ARL	1413-ARL	1414-ARL
Black	Short	1401-ARS	1402-ARS	1403-ARS	1404-ARS
	Regular	1401-ARR	1402-ARR	1403-ARR	1404-ARR
	Long	1401-ARL	1402-ARL	1403-ARL	1404-ARL



ADVANCE Armsleeve in Black

ADVANCE Armsleeve in Beige

Material composition
50% Nylon, 28% Modal, 22% Elastane (content may vary slightly by compression class)
Latex free



Scan here for this product's measuring form

ADVANCE Armsleeve

Offers breathable, moisture-wicking comfort with the option of a Sensinnov grip top-band to keep the sleeve securely in place.

Compression class

14-18 mmHg	20-25 mmHg	30-40 mmHg
CLASS 1	CLASS 2	CLASS 3

Product benefits

Moisture-wicking	Grip top band	Breathable

Suitable for

- Patients with arm lymphoedema
- Post-mastectomy care

EMIS search tips

Product name
ADVANCE Armsleeve

SIG ADV
EMIS TIP

Compression class (mmHg level)

Class 1 14 EMIS TIP Class 2 20 EMIS TIP Class 3 30 EMIS TIP

Sizing chart

All measurements are shown in cm

Armsleeve


	XS‡		S		M		L		XL	
Palm (a)	16-20		17-21		18-22		19-23		20-24	
Wrist (c)	13-14		14-16		16-18		18-20		20-23	
Upper mid arm (f)	STD	PLUS*	STD	PLUS*	STD	PLUS*	STD	PLUS*	STD	PLUS*
	24-32	32-38	26-34	34-40	29-37	37-43	32-40	40-46	34-42	42-48
EMIS TIP	SMA	SMA PLU	SMA	SMA PLU	MED	MED PLU	LAR	LAR PLU	LAR	LAR PLU

‡XS size not available in Compression Class 3 *Plus sizes only available with grip top

Length

	Normal	Long
Arm (C-G)	≤ 43	≥ 43
EMIS TIP	NOR	LON

Armsleeve



Styles and colours

- With hand piece HAN
- Without hand piece

	Class 1		Class 2		Class 3
	Black	Natural	Black	Natural	Natural
Armsleeve without grip top	● ○	○ ○	● ○	○ ○	-
Armsleeve with grip top	● ○	○ ○	● ○	○ ○	○ ○
EMIS TIP	BLA	BEI	BLA	BEI	NAT




Scan here for this product's ordering codes

ADVANCE Armsleeve



Product codes

Arm without grip top




XS	STD	Normal	Class 1				Class 2			
			Black	Black	Natural	Natural	Black	Black	Natural	Natural
			●	○	●	○	●	○	●	○
XS	STD	Normal	54714-A	54674-A	54653-A	54613-A	54724-A	54684-A	54663-A	54623-A
		Long	54719-A	54679-A	54658-A	54618-A	54729-A	54689-A	54668-A	54628-A
S	STD	Normal	54715-A	54675-A	54654-A	54614-A	54725-A	54685-A	54664-A	54624-A
		Long	54720-A	54680-A	54659-A	54619-A	54730-A	54690-A	54669-A	54629-A
M	STD	Normal	54716-A	54676-A	54655-A	54615-A	54726-A	54686-A	54665-A	54625-A
		Long	54721-A	54681-A	54660-A	54620-A	54731-A	54691-A	54670-A	54630-A
L	STD	Normal	54717-A	54677-A	54656-A	54616-A	54727-A	54687-A	54666-A	54626-A
		Long	54722-A	54682-A	54661-A	54621-A	54732-A	54692-A	54671-A	54631-A
XL	STD	Normal	54718-A	54678-A	54657-A	54617-A	54728-A	54688-A	54667-A	54627-A
		Long	54723-A	54683-A	54662-A	54622-A	54733-A	54693-A	54672-A	54632-A

● With hand piece ○ Without hand piece

Product codes

Arm with grip top



XS	STD	Normal	Class 1				Class 2				Class 3	
			Black	Black	Natural	Natural	Black	Black	Natural	Natural	Natural	Natural
			●	○	●	○	●	○	●	○	●	○
XS	STD	Normal	53835-A	53815-A	49890-A	33499-A	53845-A	53825-A	49900-A	33512-A	-	-
		Long	53840-A	53820-A	49895-A	33504-A	53850-A	53830-A	49905-A	33517-A	-	-
XS	PLUS	Normal	53855-A	54694-A	51506-A	54633-A	53865-A	54704-A	51516-A	54643-A	-	-
		Long	53860-A	54699-A	51511-A	54638-A	53870-A	54709-A	51522-A	54648-A	-	-
S	STD	Normal	53836-A	53816-A	49891-A	33500-A	53846-A	53826-A	49901-A	33513-A	63602-A	63586-A
		Long	53841-A	53821-A	49896-A	33505-A	53851-A	53831-A	49906-A	33518-A	63606-A	63590-A
S	PLUS	Normal	53856-A	54695-A	51507-A	54634-A	53866-A	54705-A	51517-A	54644-A	63610-A	63594-A
		Long	53861-A	54700-A	51512-A	54639-A	53871-A	54710-A	51523-A	54649-A	63614-A	63598-A
M	STD	Normal	53837-A	53817-A	49892-A	33501-A	53847-A	53827-A	49902-A	33514-A	63603-A	63587-A
		Long	53842-A	53822-A	49897-A	33506-A	53852-A	53832-A	49907-A	33519-A	63607-A	63591-A
M	PLUS	Normal	53857-A	54696-A	51508-A	54635-A	53867-A	54706-A	51518-A	54645-A	63611-A	63595-A
		Long	53862-A	54701-A	51513-A	54640-A	53872-A	54711-A	51524-A	54650-A	63615-A	63599-A
L	STD	Normal	53838-A	53818-A	49893-A	33502-A	53848-A	53828-A	49903-A	33515-A	63604-A	63588-A
		Long	53843-A	53823-A	49898-A	33507-A	53853-A	53833-A	49908-A	33520-A	63608-A	63592-A
L	PLUS	Normal	53858-A	54697-A	51509-A	54636-A	53868-A	54707-A	51519-A	54646-A	63612-A	63596-A
		Long	53863-A	54702-A	51514-A	54641-A	53873-A	54712-A	51525-A	54651-A	63616-A	63600-A
XL	STD	Normal	53839-A	53819-A	49894-A	33503-A	53849-A	53829-A	49904-A	33516-A	63605-A	63589-A
		Long	53844-A	53824-A	49899-A	33508-A	53854-A	53834-A	49909-A	33522-A	63609-A	63593-A
XL	PLUS	Normal	53859-A	54698-A	51510-A	54637-A	53869-A	54708-A	51520-A	54647-A	63613-A	63597-A
		Long	53864-A	54703-A	51515-A	54642-A	53874-A	54713-A	51526-A	54652-A	63617-A	63601-A

● With hand piece ○ Without hand piece



TRADITIONAL Armsleeve with shoulder cap



TRADITIONAL Armsleeve hand piece



TRADITIONAL Armsleeve with grip top

Material composition

70% Nylon, 30% Natural rubber latex (content may vary slightly by compression class)



Scan here for this product's measuring form

TRADITIONAL Armsleeve

Soft yet durable, offering reliable support and comfort, with options for a grip top or shoulder cap for a secure fit.

Compression class

23-32 mmHg

RAL 2

Product benefits



Natural rubber latex



Durable



Containment

Suitable for

- Stage I and II lymphoedema
- Post-mastectomy care

EMIS search tips

Product name
TRADITIONAL Armsleeve

SIG TRA
EMIS TIP

Sizing chart

All measurements are shown in cm


Armsleeve

	S	M	L
Wrist (C)	13-16	16-19	19-22
Mid forearm (D)	20-24	24-28	28-32
Upper Mid Arm (F)	24-30	30-36	36-42
EMIS TIP	SMA	MED	LAR

Length

	Short	Long
Arm (L _G)	40-46	46-52
EMIS TIP	NOR	LON

Armsleeve



Styles and colours

- With hand piece
- Without hand piece

HAN
EMIS TIP

Class 2

Beige

Armsleeve with shoulder cap	ARM	● ●
Armsleeve with grip top	ARM GRI	● ●
EMIS TIP	BEI	




Scan here for this product's ordering codes

TRADITIONAL Armsleeve



Product codes


Armsleeve with shoulder cap



		Class 2	
		Beige	Beige
		○	●
S	Short	80966	80972
	Long	80967	80973
M	Short	80968	80974
	Long	80969	80975
L	Short	80970	80976
	Long	80971	80977

- With hand piece
- Without hand piece

Armsleeve without shoulder cap



		Class 2	
		Beige	Beige
		○	●
S	Short	80978	80984
	Long	80979	80985
M	Short	80980	80986
	Long	80981	80987
L	Short	80982	80988
	Long	80983	80989

- With hand piece
- Without hand piece



Over **60%** of patients **struggle** with compression garments due to pain, stiffness or reduced hand strength⁸

Donning and doffing aids

Products available

DOFF N'DONNER ROLLY	86-87
DOFF N'DONNER CONE	
MELANY	



Scan the QR code to learn how to use our donning and doffing aids

DOFF N' DONNER ROLLY

Features and benefits

- Suitable for applying garments for below knee, thigh and armsleeve.
- Unlike many other donning and doffing aids, the DOFF N'DONNER ROLLY is designed to make it easier to remove and apply the compression garment.



EMIS TIP

ROLL

DOFF N' DONNER CONE

Features and benefits

- The Cone facilitates the easy rolling of compression garments onto the DOFF N' DONNER ROLLY device, streamlining the donning process.
- Equipped with a suction cup at its base, the DOFF N'DONNER CONE securely attaches to clean, flat surfaces, providing stability during use.

EMIS TIP

DOF DON CON



Using the DOFF N'DONNER ROLLY and DOFF N'DONNER CONE in conjunction makes application simple

Features and benefits



Quick and easy donning and doffing of the compression garments.



Comfortable and convenient, even with scars or ulcerated legs.

Open skin areas and wounds should be covered with an appropriate wound dressing before applying the compression garment.



Increased independence as there is less effort required application and removal of compression garments.



Use of the DOFF N' DONNER products can prevent damage to the garments as they are handled more gently.



Optimal support of nursing staff in everyday routines with compression garments.



Optimal handling for **easy cleaning**.

MELANY

Features and benefits

- Less effort required and therefore less physical strain.
- Simplified donning and doffing increases independence.
- Reduces wear and tear of compression stockings due to gentle donning and doffing allowing them to retain their therapeutic compression for longer.
- Is easy to take apart and stow away for travel or general storage.



EMIS TIP

SIG MEL


Sizing chart

All measurements are shown in cm

Below knee

	XS	S	M	L / L Plus
Calf (cC)	25-30	30-38	38-46	46-56

Below knee



Scan the QR code for the accessories form

DONNING AND DOFFING AIDS

Product codes

MELANY

MELANY

Calf

	XS	S	M	L
Product Code	309386	309387	309388	309389



DOFF N'DONNER ROLLY

One size fits all
64655



DOFF N'DONNER CONE

One size fits all
301371



Medical terms

Chronic Venous Disease	A disease state comprised of disorders ranging from tired, achy legs, varicose veins, skin issues to venous leg ulcers (open and healed).
Deep Vein Thrombosis (DVT)	A blood clot in a deep vein, most commonly in the lower leg or thigh.
Lipoedema (Lipalgia)	Lipoedema is a long-term (chronic) condition where there's an abnormal build-up of fat cells in the arms or legs.
Lymphoedema	Lymphoedema is a long-term (chronic) condition that causes swelling in the body's tissues. It can affect any part of the body, but usually develops in the arms or legs.
Lymphorrhoea (wet legs)	Lymph fluid leaking from the skin.
Mixed Diseases	A combination of venous / arterial insufficiency.
mmHg (millimetres of Mercury)	Defined as the extra pressure generated by a column of mercury one millimetre high; the same measurement used for blood pressure.
Oedema	Oedema is fluid retention in your body. The build-up of fluid causes affected tissues to become swollen.
Post-Mastectomy	Can occur over a period of time after a mastectomy procedure.
Post-Thrombotic Syndrome	Post-thrombotic syndrome is a condition that can happen to people who have had a deep vein thrombosis (DVT) of the leg. The condition can cause chronic pain, swelling, and other symptoms in your leg.
Varicose Veins	Varicose veins are enlarged, swollen, and twisted veins, often appearing blue or dark purple. They happen when faulty valves in the veins allow blood to flow in the wrong direction or start to pool.

Product material

Aramid Kermel	This material helps to improve the strength, durability and performance of the product.
Elastane	Elastane is the generic term used to describe branded textiles such as Lycra, This type of fabric is also called spandex, and its primary attribute is its incredible elasticity.
Microfibres	An extremely fine synthetic fibre, knitted into textiles with the texture and drape of natural-fibre fabric. This fabric has excellent breathability.
Modal	Modal is a semi-synthetic fibre that is used as a cotton alternative. It is both breathable and absorbent.
Neoprene	Neoprene helps heat retention, which in turn increase the blood circulation to promote proper healing.
Polyamide	A general term that is used to refer to a variety of different natural or artificial fabrics that are made from strings of polyamide monomers. The most famous form is nylon.
Polyamide Aquarius	The fabric gives an extra soft feel and a natural second-skin effect.
Polyurethane foam	Flexible polyurethane foam help in providing comfort all day.
Polyethylene	Light, versatile synthetic resin made from the polymerisation of ethylene.
Rayon	A manufactured fibre which can imitate the feel and texture of some natural fibres.
Sensinnov	A patented Sigvaris silicone top band that gently grips the skin and keeps the band in place all day.
Thermoregulating	Is a process that allows your body to maintain its core internal temperature. All thermoregulation mechanisms are designed to return your body to homeostasis. This is a state of equilibrium.

SIGVARIS GROUP

Deciding the **right product** for the **right condition** at the **right time.**

CARE AND COMPRESS

2025

DECIDING THE RIGHT PRODUCT, FOR THE RIGHT CONDITION AT THE RIGHT TIME.



Scan the QR code to find out more around Care and Compress

With over 160 years of expertise within medical compression wear, SIGVARIS GROUP is dedicated to help people feel their best. Every day



References

1. 'The treatment of venous legs ulcers with a specially designed compression stocking kit compression with bandaging.' F. Mariani et al. *Phlebologie* 2008; 37: 191-197
2. Daniel O, Mesharck G, Otor C. Effects of stress and burnout among NHS adult nurses in the UK: a systematic literature review. *J Med Public Health*. 2024;5(1):1098. Available from: <https://www.medtextpublications.com/open-access/effects-of-stress-and-burnout-among-nhs-adult-nurses-in-1592.pdf>
3. Royal College of Nursing. *Employment Survey 2019*. London: RCN; 2019. Available from: [https://www.rcn.org.uk/-/media/royal-college-of-nursing/documents/publications/2019/november/007-927.pdf​::contentReference\[oaicite:0\]\[index=0\]](https://www.rcn.org.uk/-/media/royal-college-of-nursing/documents/publications/2019/november/007-927.pdf​::contentReference[oaicite:0][index=0])
4. Mościcka P, Szewczyk MT, Cwajda-Białasik J, Jawień A. The role of compression therapy in the treatment of venous leg ulcers. *Adv Clin Exp Med*. 2019;28(6):847–852. doi:10.17219/acem/78768.
5. Damstra RJ, Partsch H. Prospective, randomized, controlled trial comparing the effectiveness of adjustable compression Velcro wraps versus inelastic multicomponent compression bandages in the initial treatment of leg lymphedema. *J Vasc Surg Venous Lymphat Disord*. 2013;1(1):13–19. doi:10.1016/j.jvsv.2012.05.001.
6. Reissbauer A, Tsafrailis E, Kornappel S, Liebl ME. A randomized controlled non-inferiority study of adjustable compression wraps compared with inelastic multilayer bandages for the treatment of lower limb lymphedema. *J Vasc Surg Venous Lymphat Disord*. 2025;13(2):123–130. doi:10.1016/j.jvsv.2024.10.005.
7. Quereda A, et al. *District Nursing Today: The Views of Team Leaders Revealed in QNI Report*. Queen's Nursing Institute; 2024.
8. Rabe, E., et al. (2020). *British Journal of Nursing*, 29(8), S30–S36.

SIGVARIS
GROUP

SIGVARIS Britain Ltd

1 Imperial Court, Magellan Close,
Andover, Hampshire, SP10 5NT

CUSTOMER SERVICE TEL: 01264 326 666

CUSTOMER SERVICE FAX: 01264 369 950

WEB: sigvaris.com



Find us on social media
SIGVARIS GROUP Britain